Telephone: LANGHAM 0157 & 0158 Telegraphic Address: LITHOTRITE, PHONE, LONDON

JOHN WEISS & SON, LTD.

Surgical Instrument Makers,

MAPPIN HOUSE WINSLEY ST., OXFORD ST., LONDON, W.1 ESTABLISHED 1787

HAND OPERATING LAMP



Complete in carrying case with transformer for A.C. Mains supply—£23 10s.

CLIFFORD BROWN L. 45 WIGMORE STREET LONDON, W.I

Telephone WELBECK 4730 (Private Branch Exchange)

Manufacturing Dispensing Opticians

A High Grade Dispensing Service is available at all our Branches

EAST SHEEN -	-	210 Upper Richmond Road, S.W.14
GLASGOW -	-	141 Bath Street, C.2
KINGSTON -	-	14 High Street
MANOR PARK	-	652 Romford Road, E.12
MORDEN -	-	8 Morden Court Parade
ORPINGTON, KEN	IT.	260 High Street
PERTH	-	143 South Street
PUTNEY	-	151 Upper Richmond Road, S.W.15
REDHILL	-	113 Station Road
REIGATE	-	25 Church Street
RICHMOND -	-	- 7 Lichfield Terrace, Sheen Road
S. WIMBLEDON	-	70 The Broadway, S.W.19
SUNDERLAND	- '	204 High Street West
TWICKENHAM	-	24 York Street
VICTORIA -	-	52 Sloane Square, S.W.I
WANDSWORTH	-	67 East Hill, S.W.18
WATERLOO -	-	298 Waterloo Road, S.E.I
WELLING, KENT	-	35 Bellegrove Road
WIMBLEDON -	-	84 High Street, S.W.19
WORTHING -	-	- Revenue Buildings, Chapel Road



L.V.A.




SCHEME

This scheme is a planned and comprehensive solution to a hitherto insoluble problem. The three sections of the scheme are:

I. ASSESSMENT OF MAGNIFICATION - A Quantitative Method.

The new KEELER "A" Series Test-Type notation, laying special emphasis on the region beyond 6/36 provides a uniform progression of type sizes. Thus, for the first time, the patient's reading acuity can be accurately determined, and — more important still — the precise magnification needed to raise the acuity to any desired or practicable level can be read off from a specially-devised table.



2. PRACTICAL VERIFICATION.

Because magnification can now be predicted without recourse to a magnifier, it follows that verification of the finding becomes a distinct and valuable procedure. The Verification Set LVA I–I shown here provides magnifiers from $2\times$ to $8\times$ all interchangeable in a trial frame, plus illuminated compound aspheric magnifiers to extend the range from $8\times$ to $20\times$.

3. DISPENSING.

A comprehensive set LVA 3-I enables a wide variety of patient's aids to be made up easily from interchangeable, standardised components. Any desired type of appliance, incorporating the magnification prescribed by the oculist, may be assembled for "on-the-spot" trial by the patient, before finally ordering.

Examples of such appliances are illustrated on facing page.



C. DAVIS KEELER LTD.

or helping the partially sighted





CRESCENT MAGNIFIER 2, 3, 4, 5 or 6x

vides distance vision while maining full width of field for reading.

TELESCOPIC SPECTACLE 1.75 or 2.5x


A monocular or binocular unit for distance vision. The decentred telescopes ensure adequate lower peripheral field.

SPECTACLE MAGNIFIER 6x or 8x

ding matter is kept in constant p focus by locating it against the sparent plastic distance-piece.

Displayed components, below show how an $R \times .$ lens may be mounted behind the aspheric magnifier.






39, WIGMORE ST., LONDON, W.I

WELBECK 8512



Melson Wingate Ltd

7, WIGMORE STREET, LONDON, W.1.



Patent applied for

THE TEMPLE GRIP TRIAL FRAME

Designed to fit all faces.


Allows for asymmetrical eyes.

Only one screw adjustment to each eye.

THE TEMPLE GRIP FRAME

A Spectacle that does not rest on the nose, yet retains the appearance of a normal frame.

Especially useful in cases of patients who cannot tolerate any pressure on the bridge of the nose.



ACCESSORIES FOR USE WITH THE SLIT-LAMP



DESIGNED AND CALIBRATED
IN THE LABORATORIES OF
THE INSTITUTE OF
OPHTHALMOLOGY.

CORNEAL PACHOMETER

FOR DETERMINATION OF CORNEAL THICKNESS

A DEVELOPMENT OF THE INSTRUMENT DESCRIBED BY MAURICE AND GIARDINI (BRIT. J. OPHTH. 1951, 35, 169—177).

THE SLIT BEAM IS DIVIDED INTO TWO RELATIVELY MOBI PARTS BY A ROTATABLE PERSPEX BLOCK. THE TWO SECTIONS OF THE CORNEAL IMAGE ARE MOVED UNITHE ENDOTHELIUM OF ONE COINCIDES WITH TI EPITHELIUM OF THE OTHER WHICH REMAINING ST TIONARY IS COLOURED BLUE BY THE INCORPORATIC OF A THIN FILTER, LOW LIGHT INTENSITY RESTRIC THE BLUE REFLEX TO THE EPITHELIUM.

THIS INSTRUMENT IS PRODUCED ON BEHALF OF STERCKS MARTIN LTD.
BY A WELL-KNOWN BRITISH MANUFACTURER.

PRICE £18 18s.

SECOND OBSERVER EYEPIECE

MANUFACTURED BY REICHERT, THE AUSTRIAN INSTRUMENT FIRM, THIS EYEPIECE IS LISTED BY HAAG-STREIT FOR USE WITH THE GOLDMANN SLIT-LAMP.

PRICE £28 7s.





GONIOSCOPY CONTACT LEN

BY GOLDMANN

AVAILABLE WITH WIDE OR MINIMAL SCLERAL FLANG (THE INSTRUMENT IS FULLY DESCRIBED IN DESCRIBED IN DESCRIBED OF DESCRIPTION OF DES

PRICE £10 10s.

AS OFFICIALLY APPOINTED AGENTS FOR HAAG-STREIT



GOLDMANN SLIT - LAMP

> AND BIOMICROSCOPE

3ABLY THE MOST POPULAR SLIT-LAMP IN USE TODAY THE BEAUTIFULLY FINISHED G-STREIT INSTRUMENT ACHIEVES BY DIRECTNESS OF APPROACH TO ALL ASPECTS BIOMICROSCOPY WITH EXTREME SIMPLICITY OF DESIGN AND STURDINESS OF ISTRUCTION THE DUAL IDEAL OF FUNCTIONAL FULFILMENT AND ECONOMY.

AVAILABLE FROM STOCK £275.

IER DETAILS OR


NSTRATION OF THE

AMP AND ACCESSORIES

PLICATION.



NELCOME ENQUIRIES REGARDING THEIR PRODUCTS





CLEMENT CLARKE

SERVICE THROUGH FORTY BRANCHES

Presenting

A NEW PAD BRIDGE FRAME

specially designed for children




Fitted with safety lenses, this frame is eminently suitable for patients old enough to be concerned with their appearance, and young enough for their parents to be concerned with the safety of the spectacles provided.

The main structure of the frame is of stainless steel, the rims and sides being covered in plastic which is available in a variety of attractive colours. The sides are curled, rigid or flexible, as required.

CLEMENT CLARKE

SERVICE THROUGH FORTY BRANCHES



BORALINE

Sole Distributors :
COATES & COOPER LTD.
WEST DRAYTON, MIDDLESEX

decongestive mildly antiseptic ideal in conjunctival irritation

OCULAR PROSTHETICS

MR. STEER WARDMAN

Artificial Eye Specialist

Brooks & Wardman

Dispensing Opticians

Artificial Eye Makers

24 REGENT STREET **NOTTINGHAM**

Telephone: 46309 2 Lines


Also attends by appointment at the following addresses:

C. DAVIS KEELER, LTD. 55. Cornwall Street. BIRMINGHAM.

Phone: Central 5422-3

J. W. WARDMAN LTD. 76, Skinnergate, DARLINGTON. Phone: 3611

BROOKS & WARDMAN LTD. 28. Watson Road. WORKSOP, NOTTS. Phone: 3204



My work covers all recent advances in ocular prostheses, magnetic artificial eyes, and facial prostheses.



DISPENSING OPTICIANS

Members of the Association of Dispensing Opticians

ANTERIOR CHAMBER IMPLANTS

We are manufacturing acrylic implants for the anterior chamber.

Our process makes it possible for wide variations in design to be provided.

We have had the advantage of discussion with eminent authorities on this subject, and patterns already manufactured are based upon their recommendations.

We shall be very pleased to send you full details, including those of a complete sterilisation service.

100 NEW BOND STREET, LONDON, W.I

Telephones: GROsvenor 5081/2/3



BY APPOINTMENT OPTICIANS

TO HIS ROYAL HIGHNESS


TO H. M. QUEEN ELIZABETH

TO THE LATE

THEODORE HAMBLIN

DISPENSING OPTICIANS
15, WIGMORE STREET, LONDON, W.1.

HAMBLIN PORTABLE EYE MAGNET



Hitherto, high power in a magnet has been possible only at the expense of mobility. This highly-efficient electro-magnet, though small in dimensions and light in weight, is capable of exerting a tractive force comparable with well-known giant types. Its shape and balance make possible a very precise approach to the foreign body and enable the operator to make the maximum use of the very high power developed.


The length (without applicator) is $5\frac{1}{4}$ inches, the maximum diameter 4 inches and the weight 6 lbs; the loading is 500 watts.

With such a loading in a magnet of these dimensions, considerable heat is generated. This is reduced to a minimum by the employment of a fin radiator.

Operating-theatre regulations require flame-proof switchgear to be used with this type of apparatus.

Where Alternating Current only is available, a rectifier is necessary. This can be supplied with suitable footswitch.

For Direct Current a control board with footswitch is also available.



PARKE, DAVIS & Company Ltd. (Inc. U.S.A.)

*Trade mark

HOUNSLOW · MIDDLESEX · Telephone : Hounslow 2361

336A xvii







BY APPOINTMENT OPTICIANS

TO HIS ROYAL HIGHNESS

TO H M. QUEEN ELIZABETH THE QUEEN MOTHER TO THE LATE

THEODORE HAMBLIN

DISPENSING OPTICIANS

HAMBLIN PURVIS STREAK RETINOSCOPE

The Purvis Streak Retinoscope made by Hamblin is of special value because it produces a reflex less confusing than that from a round spot of light, and because it enables the axis of astigmatism to be located with accuracy.

The strength of the light is controlled by a rehostat incorporated in the switch.

The width of the beam or streak of light is controlled by the sliding sleeve. The sleeve is pushed up to make the beam convergent and down to make it divergent. The divergent beam corresponds to the use of a plane mirror, and the convergent beam to a concave mirror.


For plane mirror Retinoscopy, it is only necessary to move the sleeve downwards, past the point at which the light is in sharp focus at one metre, until the beam is from $\frac{1}{4}$ to $\frac{1}{2}$ inch in width.

Alterations in the axis of the streak are governed by rotating the lower sleeve around the stem of the instrument by a lateral rolling movement of the forefinger.

THEODORE HAMBLIN [12] DISPENSING OPTICIANS 15 WIGMORE STREET, I ONDON,W.1.



For a quick, easy and accurate VISUAL FIELD testing



- BOWL made of transparent plastic. White screen painted on the inside; meridians and parallels engraved on the outside. Illumination uniform and adjustable.
- ▲ Ever-clear TEST OBJECT projected from the outside, the position of which is easy to mark.
- During the testing, the limits of the visual field are designed straight on the bowl. After examination, they are recorded on special tracing paper CHARTS. Scotomas can be recorded apart, directly traced on the bowl, with the accuracy of a CAMPIMETER.



Advertised and Introduced ONLY to the Medical Profession.


BENGUE'S

GUTTILIN

AUTOMATIC DROPPER TUBES

(BENGUE PATENT)

- ★ Sterile solutions for Ophthalmic, Aural and Nasal medication.
- ★ The Guttilin tube prevents contamination of solution.
- ★ Precision and ease of application of drops by either Practitioner or patient.
- ★ Stability and Therapeutic activity of solutions are maintained indefinitely.
- ★ Large range of solutions available.
- ★ Prices very reasonable. Economical to use, no wastage owing to precision.



List of available solutions and prices on request.

All Guttilin tubes contain 12 c.c. of solution except those of Cocaine solutions (10 c.c.) and Chloramphenicol (5 c.c.).

BENGLE & CO. LID. Mindales

MOUNT PILASANT ALPERTON WEMBLES, MEDOLISEN



Demonstration of Artificial Eye-Making given at INTERNATIONAL OPHTHALMOLOGICAL CONGRESS EXHIBITION



Artificial Eyes (Glass or Plastic) Embodying All Modern Improvements.

Specially Made to Surgeons' Instructions and Patients' Requirements

Appointments may be made at:-

2 BLOOMSBURY PLACE, LONDON, W.C.I.

(Nearest Underground Station, Holborn)
Telephone: Mus. 0908

372A

Eye Drops \ Eye Ointment

CORTISYL

1% 3ml. Bottle

1% 3G. Tube

HydroCORTISYL HYDROCORTISONE

1% 3 ml. Bottle

2.5% 3G. Tube



The Full Range

DeCORTISYL PREDNISONE

Tablets

PreCORTISYL PREDNISOLONE

Tablets

CORTISYL CORTISONE

Tablets and Injection

HydroCORTISYL HYDROCORTISONE

Skin Lotion and Ointment Tablets and Local Injection

Intravenous Injection



London, N.W.10 LADbroke 3608

LS/20

Glaxo Cortisone Products SYSTEMIC LOCAL 22 OPHTHALMIC 23

For your information . . . post free

Glaxo's "Practical Guide to Cortisone Treatment of Eye Diseases"

At your request a copy of this book will be sent to you post free. You will find it a very useful guide to all aspects of treatment with Glaxo cortisone preparations.

Glaxo Ophthalmic Preparations

Cortelan Eye Drops 1%
Cortelan Eye Ointment 1%
Ef-Cortelan Eye Drops 1%
Ef-Cortelan Eye Ointment 2.5%

Ef-Cortelan with Neomycin Eye Drops

Ef-Cortelan with Neomycin Eye Ointment

EF-CORTELAN

CORTELAN


(hydrocortisone acetate or alcohol Glaxo)

(cortisone acetate Glaxo)

GLAXO LABORATORIES LTD., GREENFORD, MIDDLESEX. BYRON 3434

Makers of ALL that's needed for Cortcosteroid treatment.

380A XXIII



The latest

RODENSTOCK EQUIPMENT

REFRACTOMETER REFRACTIONOMETER COMBI UNITS, Etc.

Send for particulars from

STEARMAN OPTICAL CO. LTD.

7, FITZROY SQUARE, LONDON, W.I Telephone: EUSton 3396

Ocular PROSTHESIS

NUMBER I O



Another of the many children successfully fitted by the special moulding methods of—

Right Eye enucleated at the age of 3 after injury. For three years badly fitting stock eyes caused shrinkage of the socket and deformed the lids. The photograph shows complete restoration of socket and lids, the result of four years of progressive work

Mr. ALFRED LESSER

A master of world-wide repute in Ocular and Facial Prosthesis.

12a, ETON AVENUE, LONDON, N.W.3. Telephone, SWIss Cottage 4949

	Y (\sim		0
A	,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	XV.	\boldsymbol{A}	j

DUBLIN

BRUSSELS

PARIS
140 Av. Victor Hugo, Paris, XVIs

41 Grafton Street, Dublin Tel: : 73741

137 Avenue Voltaire, Brussels Tel: 152278

Tel: Kleber 0060

.....


XXV

24

384A

HAWES READING SCREEN

With Adjustable Aperture



A SIMPLE DEVICE WHICH HAS PROVED VERY HELPFUL TO ELDERLY PATIENTS WITH INCIPIENT SENILE CATARACT OR MACULA CHANGES. EXTRANEOUS LIGHT IS ABSORBED BY THE SCREEN. THE APERTURE IS ADJUSTED TO A SUITABLE WIDTH LEAVING A FEW LINES OF PRINT CLEARLY DISCERNIBLE.

ALFRED HAWES & SON, LTD.

79, LEADENHALL STREET, LONDON, E.C.3 25, New Cavendish Street, London, W.1


and at

CHISLEHURST ELTHAM LEWISHAM STRATFORD

CROYDON GODALMING LEYTON WHITECHAPEL

DAGENHAM GUILDFORD ROMFORD WALTHAMSTOW WOKING WALTON-ON-THAMES

DORKING **ILFORD** ST. PAUL'S CRAY



xxvii



BY APPOINTMENT OPTICIANS

TO HIS ROYAL HIGHNESS


TO H. M. QUEEN ELIZABETH

TO THE LATE

THEODORE HAMBLIN 🖭

DISPENSING OPTICIANS

HAMBLIN UNIVERSAL SLIT-LAMP



15 WIGMORE STREET, LONDON, W.1.

=

LINWOOD STRONG

DISPENSING OPTICIANS

in

SUTTON
20 HIGH STREET

EPSOM
45 UPPER HIGH STREET

&

WORTHING
96a MONTAGUE STREET

UNIVERSITY OF LONDON

INSTITUTE OF OPHTHALMOLOGY, JUDD STREET, LONDON, W.C.I

Associated with MOORFIELDS EYE HOSPITAL

Qualified Medical Practitioners may enter the practice of the Moorfields Eye Hospital at any time, and are, on certain conditions, eligible for appointment as Clinical Assistants.

Systematic courses suitable for Diploma and Degree examinations begin in OCTOBER and MARCH of each year.

For further particulars apply to the Dean.

Examination of Pathological Ophthalmic Specimens.

The Department of Pathology of this Institute will examine and report on certain specimens sent by Ophthalmic Surgeons.

For particulars apply to the Director of the Department of Pathology.

ROYAL EYE HOSPITAL

ST. GEORGE'S CIRCUS, S.E.1


Telephone: WATERLOO 4477

The clinical practice of the Hospital is open to postgraduate students. There are 176,000 total outpatient attendances and over 2,000 admissions annually.

Courses for the Fellowship in Ophthalmology are held as required. Further information can be obtained from the Secretary.

Canada ophthalmic practice very large for sale to two ophthalmic surgeons. Introduction 6 to 12 months. Purchase office equipment only requirement. Total cost deductible from income tax. Equals six months net. Up to one half accepted in blocked Sterling (5,000 pounds). Balance by low interest five year loan. Reply to Box No. 50 British Journal of Ophthalmology, B.M.A. House, Tavistock Square, London, W.C.I.

Fashion Frames



A further example of the spectacle-maker's craft, combined with the jeweller's. Hand-made frames, trimmed hand-set rhinestones and pearls.

Earrings to match.

Made and Designed by

A. E. NEWBOLD, LTD.

11b, Hay Hill, Berkeley Square, London, W.1

HYDE PARK 0945/6

Branch: 1 Broadway Buildings BROMLEY, Kent

Branch: 64 High Street BECKENHAM, Kent

GEORGE ODELL LTD.

DISPENSING OPTICIANS

A. H. ODELL, Managing Director. M. ODELL, Secretary.

Members of the Guild of British Dispensing Opticians

Makers of spectacles of a fine precision and of a standard worthy of the best work of British craftsmen—a very important but fast dwindling section of the optical industry.

Frames are made and mounted to the individual measurements of each patient and both frame and lenses periodically checked during the course of production.

The finished spectacles are finally checked by a qualified dispensing optician and whenever possible, personally fitted on every patient.

31, KING STREET, ST. JAMES'S, S.W.I Whitehall 9987 36, WELBECK STREET, W.I

Welbeck 5320

Faithfully yours

PANKHURST & BARRETT LTD.

Dispensing Opticians

27, NEW CAVENDISH STREET, LONDON, W.I 250, IMPERIAL DRIVE, RAYNERS LANE, HARROW

ROMMEL HILDRETH CAUTERY, as used by JOHN FOSTER, F.R.C.S.



Supplied complete with battery in handle.

DOWN BROS. and MAYER & PHELPS Ltd.

CHURCH PATH, CHURCH ROAD, MITCHAM, SURREY 32-34, NEW CAVENDISH STREET, LONDON, W.I

and at

70, GRENVILLE STREET, TORONTO, CANADA



Blakiston Division
McGraw-Hill Publishing Company Ltd
McGraw-Hill House London EC4

ATLAS OF EYE SURGERY

R. Townley Paton, MD, FACS and Herbert M. Katzin, MD, FACS Illustrated by Daisy Stillwell

Designed for resident surgeons in departments of ophthalmology or for the general surgeon who occasionally has to perform eye operations, the approach in this new atlas is almost entirely visual. All the 464 fine and detailed drawings of surgical procedure are drawn from the point of view of the surgeon and clearly indicate intended procedures. The text is brief, and merely serves to summarize information which cannot be illustrated.

Published

240 pages

illustrated

£5 12s 6d

CURRY & PAXTON LTD.

THE NEW SLIT LAMP

by Carl Zeiss of Jena

Incorporating variable magnification

We shall be demonstrating this instrument on our stand together with other Carl Zeiss

Ophthalmic Equipment at the Oxford Congress on the 1st, 2nd and 3rd of July

CURRY & PAXTON LTD

195-199, GREAT PORTLAND STREET and 22, WIGMORE STREET, LONDON, W. I

(Head Office) (Instrument Dept.)

PRINCIPAL BRANCHES

LIVERPOOL SOUTHPORT RICHMOND (SURREY) NORWICH
GLASGOW BATH LOWESTOFT

BRISTOL NOTTINGHAM DERBY LINCOLN



Our Standard Trial Case now contains our major range of small aperture trial lenses fitted in a solidly made panelled, lift-off lid case with removable tray. The lid is secured by 2 strong spring clasps and a lock and key is provided. Can be supplied in natural mahogany or oak.

The lenses are well spaced for easy removal. The indices are clearly engraved on white plastic with black figures. Outside measurements $21\frac{1}{2}'' \times 15\frac{3}{4}'' \times 2\frac{3}{4}''$.

9 Accessories: Maddox Rod. Red & Green Glass. 2 Cross Lines. Occluder Disc. Pin Hole. Stenopaic Slit. Frosted Glass.

TOTAL Number of Lenses - 225

CURRY & PAXTON LTD.

195-199, GREAT PORTLAND STREET and 22, WIGMORE STREET, LONDON, W.

(Head Office)

(Instrument Dept.)

Principal Branches:

L[†]VERPOOL GLASGOW

SOUTHPORT

Ватн

RICHMOND (SURREY) Norwich Lowestoft


BRISTOL

NOTTINGHAM

DERBY

LINCOLN

Clement Clarke Present



for photography of the eye in black-and-white or colour

Product of the distinguished French Company, O.P.L., this superb Retinal Camera is now obtainable from Clement Clarke. It is impossible to give full details of the instrument here, but the following is a summary of its benefits.

It is suitable for use on 110 to 250 volts without special installation \bigstar Fitted with an interchangeable electrical unit \bigstar Ease of focusing during the examination due to the two reduction levers, without reversing the movement \bigstar Easy focusing during the examination with constant observation of the patient due to the higher position of the operator and the good side clearance \bigstar Reliability and complete synchronisation of the photographic mechanism \bigstar Advantages due to the use of flash and 35 mm film \bigstar Correction of all spherical and astigmatic ametropy \bigstar

A comprehensive booklet on the O. P. L. Retinal Camera will be sent on request

CHEMIENT CLARKE

INSTRUMENT DEPARTMENT: 63 WIGMORE STREET, LONDON W.1

AUTHORISED AGENTS IN BRITAIN FOR O.P.L. OF FRANCE

Clement Clarke Present



GOLDMANN APLANATION TONOMETER

The APLANATION TONOMETER, manufactured by HAAG-STREIT, is the result of several years of intensive study of tonometry in Professor Goldmann's clinic at Berne. The instrument is used in conjunction with the Goldmann Slit-lamp as is shown in the illustration. We have imported the first Aplanation Tonometer into this country and are now in a position to invite enquiries.

* OTHER HAAG-STREIT APPARATUS which can be supplied from Stock includes the Slit Lamp, Projection Perimeter, Ophthalmometer and Adaptometer. Servicing facilities (on site or at our Works) are available for all Haag-Streit instruments.

INSTRUMENT DEPARTMENT: 63 WIGMORE STREET, LONDON, W.1.

AUTHORISED AGENTS IN GREAT BRITAIN FOR HAAG-STREIT LTD



APPOINTMENT OPTICIANS

DISPENSING OPTICIANS 15, WIGMORE STREET, LONDON, W.1

SOLE AGENTS IN THIS COUNTRY FOR

ERNST GRIESHABER OF SCHAFFHAUSEN

The makers of the well-known FRANCESCHETTI Trephines and Needles, GIRARDET and CASTROVIEJO Keratoplasty Instruments, Needles, Knives, Keratomes, Refractors, Curettes, etc. Full detailed list on application.

Sizes of SUTURE NEEDLES illustrated on left



No. 80



No. 27 Vectis

No. 81



Curette



No. 40a Scalpel (Elschnig)




No. 13 Keratome



No. 31 Sharp




Knife Needle (Haab)



xxxix



- **SEGMENT** REFLECTION FREE
- **SEGMENT** INCONSPICUOUS
- **SEGMENT** $23 \times 16 - \frac{1}{2}$ mm.
- **BLANK** 60×58 mm. SHAPED
- TINTS CRUXITE A. & Ax.
- **AVAILABILITY BLANKS** SEMI-FINISHED FINISHED LENSES
- **BASES** PLUS 4: 6:8:10 D
- **ADDITIONS** 1.00 TO 3.00 IN 1 D STEPS & 3.50 & 4.00





Stockists will appreciate the large blank which removes the necessity for holding right and left semi-finished blanks.

BRITISH AMERICAN OPTICAL COMPANY LIMITED

39, HATTON GARDEN, LONDON, E.C.I.

Factories: WATFORD (Herts) KIDWELLY (S. Wales)

S. P. HOLLOWAY LTD.

PROMPT DISPENSING

WORKSHOP ON THE PREMISES

14a. HENRIETTA PLACE

(Corner of Welbeck Street)

LONDON · W.1

Telephone WELbeck 2703

DEVELOPMENTAL ABNORMALITIES OF THE EYE

By IDA MANN C.B.E., M.A., D.Sc., M.B., B.S., F.R.C.S., F.R.A.C.S.

SECOND EDITION, REVISED AND ENLARGED

March, 1957. 419 pages. 284 illustrations, 68 in full colour. 90s. net (postage, 1s. 6d. inland, 1s. overseas)

BRITISH MEDICAL ASSOCIATION

B.M.A. House, Tavistock Square, London, W.C.1
PHILADELPHIA - J. B. LIPPINCOTT COMPANY

[&]quot;A work which will be found to be indispensable to all practising ophthalmologists and sheds lustre on British ophthalmology."—BRITISH JOURNAL OF OPHTHALMOLOGY.

[&]quot;A beautiful book, well written, exquisitely illustrated. This book is a classic. A copy should be owned by every ophthalmologist. Though primarily for a reference, it makes exceedingly interesting reading from cover to cover."—AMERICAN JOURNAL OF OPHTHALMOLOGY.



GEORGE SPILLER,

Lta.

41, WIGMORE ST., W.I

Tele.: WEL. 7471-2

THE FIRST OPTICIAN IN WIGMORE STREET. ESTABLISHED 1873