


# JOHN WEISS & SON, LTD.


17 WIGMORE STREET,  
LONDON, W.1

Telephone: MUSeum 7733


## INSTRUMENTS FOR CORNEAL SUTURING



B951b




B874



B1053

- B951b.** Suture Forceps, corneal, St. Martin's, *stainless steel*.
- B874.** Forceps, corneal, Barraquer's *do.*
- B1053.** Scissors, miniature, for removal of corneal sutures, etc.



B607



B608



B609

- B607.** Needles corneal, Jameson Evans', 8, 10 and 14 mm.
- B608.** Ditto *do.* Kalt's, 12, 14 and 18 mm.
- B609.** Ditto *do.* Vogt's, 7, 10 and 14 mm.

You are cordially invited to visit our new Showrooms at 17, Wigmore Street, W.1, and inspect a comprehensive range of Ophthalmic Instruments and Apparatus.

**CLIFFORD BROWN <sup>D.</sup>**  
**45 WIGMORE STREET**  
**LONDON, W.1**

*Telephone: WEL 4730*

(Private Branch Exchange)

***Manufacturing Dispensing Opticians***

*A High Grade Dispensing Service is available  
 at all our Branches*

EAST SHEEN	- -	210 Upper Richmond Road, S.W.14
GLASGOW	- - - -	141 Bath Street, C.2
KINGSTON	- - - -	14 High Street
MANOR PARK	- - -	652 Romford Road, E.12
MORDEN	- - - -	8 Morden Court Parade
ORPINGTON, KENT	- - - -	260 High Street
PERTH	- - - -	143 South Street
PUTNEY	- - -	151 Upper Richmond Road, S.W.15
REDHILL	- - - -	113 Station Road
REIGATE	- - - -	25 Church Street
RICHMOND	- - -	7 Lichfield Terrace, Sheen Road
S. WIMBLEDON	- - -	70 The Broadway, S.W.19
SUNDERLAND	- - -	204 High Street West
TWICKENHAM	- - - -	24 York Street
VICTORIA	- - - -	52 Sloane Square, S.W.1
WANDSWORTH	- - - -	67 East Hill, S.W.18
WATERLOO	- - - -	298 Waterloo Road, S.E.1
WELLING, KENT	- - - -	35 Bellegrove Road
WIMBLEDON	- - - -	84 High Street, S.W.19
WORTHING	- -	Revenue Buildings, Chapel Road
WEST WORTHING	- -	6 Wallace Parade, Goring Road




BY APPOINTMENT OPTICIANS TO HER LATE MAJESTY QUEEN MARY

C. W. DIXEY & SON LTD.


19 WIGMORE STREET · LONDON W.1.

TELEPHONE: LANGHAM 1713


### Ophthalmic Surgical Instruments and Apparatus



D 1203



D 711



D701-2


- D 1203. Barraquer's Colibri Forceps, for corneal suturing, etc.
- D 711. Vannas' Scissors, for occlusio pupillae, anterior synechia, capsule, etc.
- D 701-2. Castroviejo's Incision Enlarging Scissors, Right & Left.

*Delivery of the above or other instruments from our comprehensive range can usually be made from stock.*

*Catalogue sent on request.*

# Melson Wingate Ltd

7, WIGMORE STREET, LONDON, W.1.



*Patent No.: 824876*

## THE TEMPLE GRIP TRIAL FRAME

Designed to fit all faces.


Allows for asymmetrical eyes.


Only one screw adjustment to each eye.

## THE TEMPLE GRIP FRAME

A Spectacle that does not rest on the nose, yet retains the appearance of a normal frame.

Especially useful in cases of patients who cannot tolerate any pressure on the bridge of the nose.





at the First European  
Congress of Ophthalmology

ATHENS APRIL 18th – 22nd

*Clement Clarke will be showing a selection of new and revised instruments at the Athens exhibition including the following:*

- 1 An entirely *new* Synoptophore with built-in automatic flashing unit and Haidinger's Brush device.
- 2 Haidinger's Brush in Free Space apparatus for the treatment of eccentric fixation.
- 3 Electronic Flash After-Image Test.
- 4 Rolson Projection Focimeter.

*Please write for further information.*

**CLEMENT CLARKE LTD.**

INSTRUMENT DIVISION, 63 WIGMORE STREET, LONDON, W.1.



## THE ELECTRONIC FLASH AFTER-IMAGE TEST

The Electronic Flash After-Image Test overcomes the major drawback encountered with all after-image instruments hitherto manufactured -- maintenance of steady fixation by the patient.

With this apparatus only momentary fixation is necessary, for the electronic flash tube provides an ample stimulus to promote a strong after-image in most patients, irrespective of the degree of suppression.

Further, the fixation point being self-luminous the test can be carried out in darkness, to avoid occurrence of peripheral fusion.

**CLEMENT CLARKE LTD.**

*INSTRUMENT DIVISION, 63 WIGMORE STREET, LONDON, W.1.*

# *Alfred Hawes*


*and Son Ltd.*

## THE HIGHEST STANDARDS OF OPTICAL DISPENSING

Our service to the Ophthalmic Medical Profession dates from the year 1840, for it was in that year that our firm was established in the City of London.

Tradition and service have, and always will play an important part of the life that is contained in the 'Square Mile', and although we have now extended our organisation to Branches in other parts of London and the Home Counties, we have been at pains to maintain at these other establishments the reputation that the City of London enjoys.

We look forward to dispensing the requirements of your patients, backed with the advantages of past experience coupled with attention to present advances.



**HEAD OFFICE:**

**79-79a, LEADENHALL STREET, LONDON, E.C.3.**

**TELEPHONE: AVENUE 4538/9**

**WEST END OFFICES: 25, NEW CAVENDISH STREET, LONDON, W.1. Telephone: WELbeck 7862**


**BRANCHES:**

<b>BOW:</b>	<i>Advance 2103</i>	<b>GUILDFORD:</b>	<i>Guildford 4911</i>	<b>ST. PAUL'S CRAY:</b>	<i>Foots Cray 4803</i>
<b>CHISLEHURST:</b>	<i>Imperial 4433</i>	<b>ILFORD:</b>	<i>Valentine 5825</i>	<b>STRATFORD:</b>	<i>Maryland 3393</i>
<b>CROYDON:</b>	<i>Croydon 5779</i>	<b>LEATHERHEAD:</b>	<i>Leatherh'd 3157</i>	<b>WALTHAMSTOW:</b>	<i>Coppermill 5206</i>
<b>DORKING:</b>	<i>Dorking 3187</i>	<b>LEWISHAM:</b>	<i>Lee Green 2245</i>	<b>WOKING:</b>	<i>Woking 944</i>
<b>ELTHAM:</b>	<i>Eltham 4325</i>	<b>LEYTON:</b>	<i>Leytonstone 2109</i>	<b>WALTON-ON-THAMES:</b>	<i>Walton 183</i>
<b>GODALMING:</b>	<i>Godalming 2088</i>	<b>ROMFORD:</b>	<i>Romford 42208</i>		



## Amsler Charts Manual

Seven charts in substantially bound book with set of 50 recording charts on tear-off pad, and full instructions for their use by Professor Amsler.



These charts, designed by Professor Marc Amsler of Zürich, make it possible to detect, localise, and analyse the visual functional disturbances (metamorphopsia, relative scotoma) which accompany the beginning and evolution of maculopathies. Many functional disturbances affecting the central and paracentral region of the field escape the usual quantitative test of central vision since it does not really bring into account the quality of the visual trouble. By means of these charts and the method and technique of examination described by Professor Amsler such qualitative measurements are easily made.

*Note.*—Spanish, French and German editions are available at same price.

**THEODORE HAMBLIN LTD**

15 Wigmore Street, London, W.1.


---

---

AN UNDENIABLE ASSET

## THE 'SECOND OBSERVER' EYEPiece

*THIS DEMONSTRATION DEVICE FOR USE WITH  
THE SLIT LAMP IS GAINING PARTICULAR  
ATTENTION IN TEACHING HOSPITALS*



The extra eyepiece provides the student with a first-hand view of procedure which is invaluable in GONIOSCOPY and APLANATION TONOMETRY instruction.

A pointer, operated by the demonstrator, is superimposed upon the field of view.

*Please write or telephone  
for further particulars*

**Stercks Martin**

DISPENSING OPTICIANS

INSTRUMENT DEPARTMENT

36, NEW CAVENDISH ST.,  
LONDON, W.1.

TEL. HUNTER 0048

---

---

SOLE AGENTS UNITED KINGDOM & EIRE FOR GAMBS OF LYON

# OCULAR PROSTHETICS

MR. STEER WARDMAN

*Artificial Eye Specialist*

*of*

***Brooks & Wardman***

*Dispensing Opticians*

*Artificial Eye Makers*

24 REGENT STREET

NOTTINGHAM


Telephone: 46309 2 Lines

*Also attends by appointment at the following addresses:*

91, Cornwall Street,  
BIRMINGHAM, 3.  
Phone: Central 0583

76, Skinnergate,  
DARLINGTON.  
Phone: 3611

28, Watson Road,  
WORKSOP, NOTTS.  
Phone: 3204




My work covers all recent advances in ocular prostheses, magnetic artificial eyes, and facial prostheses.



**OUT** with  
battery  
troubles

These latest **KEELER** electrical items are the final solution to the problem of ensuring constant bright light in battery instruments, making conventional dry cells a thing of the past.



### THE CHARGING UNITS


Single Unit (shown above)

- ★ Double wound transformer for safety.
- ★ Recharges Charge-a-Lite inserts or complete handles.
- ★ 2 models, 110v A.C. (1636E), 230/250v A.C. (1637E).

Double Unit (shown opposite)

- ★ Simultaneous charging of 2 Charge-a-Lite inserts or 2 complete handles.
- ★ 2 models—110v A.C. (1638E), 230/250v A.C. (1639E).

**IN** with **KEELER**  
RECHARGEABLE  
UNITS



**THE CHARGE-A-LITE** — a rechargeable battery insert (ref. 1631G).



- ★ Fits any Keeler battery handle and most other handles of U.2 battery size.
- ★ More than 3 hours capacity at maximum brilliance.
- ★ Controlled voltage ensuring greater bulb life.

The Keeler Charge-a-Lite and the Single Charging Unit are the only items required to convert any Ophthalmoscope and/or Retinoscope to modern rechargeable specifications.

## RECHARGEABLE HANDLES.

Made to the same high specifications as the Charge-a-Lite inserts the rechargeable handles (below and right) are complete entities and may be inserted directly into either charging unit. The handles are slim, well-balanced and have a durable satin-chrome finish.

Both handles and charging units are attractively styled for desk use, so that automatic recharging can take place between examinations.



1633 G    1632 G

MODEL 1632 G: designed for all instrument heads except the Specialist Ophthalmoscope and Streak Retinoscope.

MODEL 1633 G: especially designed for the "Specialist" and "Streak" but will accept all heads.

NOT JUST A NOVELTY . . .

a major contribution to trouble-free electrical power. An investment now ensures against battery failures in the future.




39 Wigmore Street, London W.1., England.

5536 Baltimore Av., Philadelphia 43, Pa. U.S.A.

# THE COLOURED HAPTIC A.C. IMPLANT

(For use in cases of aniridia)



The development and use of this, and one other new form of anterior chamber implant is described in *The Lancet* Jan. 9th 1960, p. 90.

We shall have samples available for inspection on our stand at the Trade Exhibition held in conjunction with the Ophthalmological Society United Kingdom Congress 1960 at Cowdray Hall from 31st March to 2nd April.

**Instrument Department,  
160 New Bond Street  
London, W.1**


MAYfair 0128

*Our complete Ophthalmological Catalogue gladly sent on request*

# Hamblin


## OPHTHALMIC LAMP

Provides small brilliant areas of light suitable for retinoscopy and ophthalmoscopy. The largest luminous portion supplies illumination for the selection of lenses and note-taking



## ANGLEPOISE BRACKET

In conjunction with the Hamblin Ophthalmic Lamp this makes a useful and inexpensive bracket, which can be placed in almost any position by a slight touch.



**THEODORE HAMBLIN LTD**


15 Wigmore Street, London, W.1.

# GLAUCOMA

*The new carbonic anhydrase  
inhibitor to reduce  
intra-ocular pressure*

## DARANIDE

inhibits the formation of  
aqueous humour  
lowers intra-ocular tension  
within one hour and maintains  
this reduction for six hours or  
more  
gives smooth control with  
relatively few side-effects  
reduces the likelihood of  
metabolic acidosis in  
comparison with other carbonic  
anhydrase inhibitors  
may often be effective when  
other carbonic anhydrase  
inhibitors have failed or have  
not been tolerated



# DARANIDE

(Dichlorphenamide)

Trade Mark

*FOR THE ORAL TREATMENT OF GLAUCOMA*

Tablets, 50 mg., in bottles of 100


The United Kingdom N.H.S. basic cost is 40s. 0d. per 100 tablets

Full details and references gladly supplied on request



**MERCK SHARP & DOHME LIMITED, HODDESDON, HERTS**





## ZONULOLYSIS in Cataract Surgery eliminates accident risks

The inherent complications of extracapsular extraction and the accidents which sometimes occur during or after intracapsular extraction are eliminated by the new technique, originated by BARRAQUER, which employs the zonulolytic action of  $\alpha$ -chymotrypsin and replaces effectively the mechanical procedures of zonulotomy currently used.

As pioneers in research and development of proteolytic enzymes, we have collaborated closely in this new technique by producing

# CHYMAR - ZON

brand of  $\alpha$ -chymotrypsin

**PRESENTATION**—each carton contains

1 vial—100 Armour units lyophilised  $\alpha$ -chymotrypsin.

1 vial—1 ml. Ringer's Solution as diluent.

1 vial—1 ml. Ringer's Solution for irrigation.

*Literature and further information on request.*




**ARMOUR PHARMACEUTICAL COMPANY LTD.**  
**EASTBOURNE, ENGLAND**

Telephone: Hampden Park 740.

Telegrams: 'Armo/lab' Eastbourne

## **OPHTHALMIC DIATHERMY**



### **Hamblin "OPHTHERM"**

*a choice of two models  
with output of identical frequency*

**Valve or  
Spark Gap**

Simple, safe and dependable, the Hamblin OPTHHERM is a diathermy designed especially for ophthalmic surgery. Proved extremely stable and efficient, it has an output adequate for every type of patient and every degree of tissue resistance. Output, shown on a milliammeter, is controlled by either hand- or foot-switch. A graduated scale enables any setting to be repeated with exactitude; a tick at one-second intervals facilitates the timing of treatment. The OPTHHERM is compact and easily portable, and its robust case completely protects the circuitry.

**THEODORE HAMBLIN LTD**

15 Wigmore Street, London, W.1.



20 High Street,  
Sutton, Surrey.

*Telephone: VIGilant 2798.*

Connaught Buildings,  
Chapel Road, Worthing.

*Telephone: Worthing 5193.*

45 Upper High Street,  
Epsom, Surrey.

*Telephone: EPSom 2300.*

57b London Road,  
Brighton.

*Telephone: Brighton 29425.*



# BORALINE

decongestive  
mildly  
antiseptic

ideal in  
conjunctival  
irritation



*Sole distributors:*

**COATES & COOPER LTD.** WEST DRAYTON, MIDDLESEX



## Zonulysin in cataract

Zonulysin is a highly purified, sterile freeze-dried Alpha-Chymotrypsin, which is particularly indicated for resolving the ciliary zonule.

Maw's Zonulysin is presented in a convenient form for this purpose, allowing easy preparation of a fresh solution for each procedure.



Zonulysin instilled to flood anterior and bathe suspensory ligament.



Suspensory ligament resolved giving easy withdrawal of crystalline lens.



Supplied in boxes of 6 ampoules, each containing 0.2 mgm. Alpha-Chymotrypsin together with 6 ampoules each containing 1.0 ml. of Saline diluent.

*Enquiries to Ethical Department*

**S. MAW SON AND SONS LIMITED**

ALDERSGATE HOUSE · BARNET · HERTFORDSHIRE

—|—  
**the Hamblin**  
**FIXATION FINDER**  
 —|—



### Parafoveal Fixation

This can readily be distinguished from foveal fixation by means of an ordinary Hamblin Lister-Morton Ophthalmoscope, the standard condenser of which has been replaced by this accessory. A blue filter, on which appears a central white circle, is incorporated in a condenser specially computed for focusing an image of the circle on the the patient's fundus.

### Procedure

When using the fixation finder, the ophthalmologist instructs the patient to look into the instrument, directly at the circle. The area used for fixation, contrasted by its white illumination against the general blue lighting and circumscribed by a dark ring, can be plainly seen ophthalmoscopically. The ophthalmoscope head, for an emmetropic patient, is set rather less than 1 cm. up the pillar of the instrument; for a myope, the head will be nearer the handle; for a hypermetrope, it will be farther up the pillar.

### a Lister-Morton Ophthalmoscope Accessory

Until complete familiarity has been gained, two minor difficulties may be experienced in the use of the fixation finder. Firstly, the patient may not immediately see the circle; in such a case, he should be told to search for it directly into the light, which should be dimmed. Secondly, the observer himself may not directly locate the circle; this trouble can be overcome by initially observing the image on an indifferent surface, such as the patient's forehead through the sight-hole of the ophthalmoscope. The hold so gained on the instrument is preserved for examining the patient.

**THEODORE HAMBLIN LTD**

15 Wigmore Street, London, W.1.

**BI-LATERAL ORBITAL  
AND  
NASAL PROSTHETIC  
RESTORATION**



Prosthetic Restoration as worn by patient



Patient without restoration



Patient wearing restoration

**CLEMENT CLARKE**

**ARTIFICIAL EYE AND PROSTHESES DEPARTMENT:—**

---

## CLEMENT CLARKE

### AESTHETIC PROSTHESES

As an adjunct to surgical practice we offer to Surgeons our modern service of plastic prosthetics, either as an interim measure during surgical correction, or for permanent rehabilitation in the cases where plastic surgical procedure is contra-indicated. In conjunction with our plastic artificial eye service some remarkable results have been obtained, as the illustrations on the opposite page clearly demonstrate.

Appointed official contractors to Ministry of Health.

### ARTIFICIAL EYES

Complete range and accessories supplied. Eyes fitted to all types of buried and semi-buried implants. Plastic Artificial Eyes specially made to an impression of the socket can be successfully fitted in the case of injury or surgical intervention. The irides are faithfully reproduced by hand pigmenting the stroma pattern and colours on a series of acrylic discs, superimposed to simulate the life and depth of the natural iris. Excellent results are obtained in fitting shell eyes to shrunken and deformed globes.



COSMETIC CONTACT PROSTHESES FITTED TO UNSIGHTED AND  
DISFIGURED EYES.

CONFORMER TREATMENT TO DISTEND CONTRACTED SOCKETS.  
OCULAR SPLINTS. IMPLANTS.

"ARRUGA" INCLUSION—TWIN POST TYPE AND CONFORMERS  
POST-OPERATIVE SHELLS—GLASS OR PLASTIC

---

## CLEMENT CLARKE

16 WIGMORE STREET, LONDON, W.1.    Tel. LANGHAM 2242

---



**CONTROL  
INFECTION  
NIGHT AND DAY**

with

**'BROLENE'**

*trade mark*

*brand*

DIBROMOPROPAMIDINE ISETHIONATE  
EYE OINTMENT

and

**PROPAMIDINE  
OPHTHALMIC SOLUTION**

These two preparations provide complementary forms of treatment for infections and superficial injuries of the eye. 'Brolene' ointment is best used at night while Propamidine Ophthalmic Solution is ideal for daytime application.

*PRESENTATIONS*

'Brolene' is supplied as a 0.15 per cent eye ointment.  
Propamidine Ophthalmic Solution: 0.1 per cent solution  
in dropper bottles.

*Further information is available on request.*

MANUFACTURED BY


**MAY & BAKER LTD**



*DISTRIBUTORS*  
PHARMACUTICAL SPECIALITIES  
(MAY & BAKER) LTD  
DAGENHAM

**M&B brand Medical Products**


MA6321 120



# THE CURPAX

VANDERBILT  
PROJECTION PERIMETER

*for quicker, easier plotting of the visual fields*



An improved design incorporating new features and  
with direct recording

## THE CURPAX VANDERBILT PROJECTION PERIMETER

*improved design permits:—*

★ More positive positioning of patient	by retractable rod and spring-loaded chinrest.
★ Use of direct reading charts	easily inserted and removed from paddle holder.
★ Direct plotting of both fields	on one chart, without removal from holder.
★ Recording at every 15 degrees	arc holds positive position at every 15°.
★ Use of standard 12 v. bulb	only the simple bayonet type bulb is changed.
★ Simple centring and focusing of light	a built-in control assures simplicity of adjustment
★ Transformer and switch built in	eliminates extra equipment on or under table.
★ Perimeter to fit any standard table	no overhang of chinrest, or need for special table.

## CURRY & PAXTON, LTD.

*Dispensing Opticians, Manufacturers of Ophthalmic Instruments*

**195-199, GREAT PORTLAND ST. and 22, WIGMORE ST., LONDON, W.1**

*Telephone: LANgham 0123.*


*Telephone: LANgham 2873.*

**Instrument Dept: 120, Albert Street, Parkway, N.W.1**

*Telephone: GULLiver 7939 and 5495*

PRINCIPAL BRANCHES—LIVERPOOL, BRISTOL, BATH, SOUTHPORT, GLASGOW, DERBY, LOWESTOFT, WARRINGTON, WESTON-SUPER-MARE, NOTTINGHAM, RICHMOND, NORWICH, LINCOLN.

## *Fashion Frames*



\* \*

*A further example of the spectacle-maker's craft,  
combined with the jeweller's. Hand-made frames,  
trimmed hand-set rhinestones and pearls.  
Earrings to match.*

\* \*

*Made and Designed by*

**A. E. NEWBOLD, LTD.**

**11b, Hay Hill, Berkeley Square, London, W.1**

*HYDE PARK 0946/6*

Branch: 1 Broadway Buildings  
BROMLEY, Kent

Branch: 64 High Street  
BECKENHAM, Kent

## ACCESSORIES FOR CORNEAL GRAFTING

### Pittar Splint

Made of thin tantalum metal, it has a central perforation for drainage and four peripheral lugs for the attachment of retaining sutures.



### Pittar Trepine Guide

The guide facilitates the centration of the trephine prior to securing a graft with a Pittar or a Stallard Splint.




### Stallard Splint

Similar to the Pittar pattern, but without perforation and made in transparent acrylic so that the entire graft may be kept under observation.



### Philps Splint

Made in acrylic, it has peripheral holes for the reception of sutures and additional perforations for drainage. An engraved cross enables the splint to be centred with ease.



### Ainslie Splint

Ainslie's Modification of the Philps Splint differs from the original only in having considerably larger drainage holes.



### Ainslie Corneal Marker

The Corneal Marker is for marking precisely on the cornea the locations for the retaining sutures. Its eight short needle points exactly correspond with the suture holes of the splint.

**Hamblin**  
15 Wigmore St.  
London W.1.

# GEORGE ODELL LTD.

## DISPENSING OPTICIANS

A. H. ODELL, *Managing Director.*  
M. ODELL, *Secretary.*

*Members of the Guild of  
British Dispensing Opticians*

Makers of spectacles of a fine precision and of a standard worthy of the best work of British craftsmen—a very important but fast-dwindling section of the optical industry.

Frames are made and mounted to the individual measurements of each patient and both frame and lenses periodically checked during the course of production.

The finished spectacles are finally checked by a qualified dispensing optician and whenever possible, personally fitted on every patient.

**31, KING STREET,  
ST. JAMES'S, S.W.1**  
Whitehall 9987

**36, WELBECK STREET  
W.1**  
Welbeck 5320

**S. P. HOLLOWAY LTD.**

**PROMPT DISPENSING**

**WORKSHOP ON THE PREMISES**


**14A, HENRIETTA PLACE**

*(Corner of Welbeck Street)*

**LONDON · W.1**

*Telephone Welbeck 2703*

### ERISIPHAKES



DOWN BROS. and MAYER & PHELPS Ltd.

CHURCH PATH, MITCHAM, SURREY  
32-34, NEW CAVENDISH STREET, LONDON, W.1  
and at  
70, GRENVILLE STREET, TORONTO, CANADA

*Faithfully yours*

**PANKHURST & BARRETT LTD.**

*Dispensing Opticians*

27, NEW CAVENDISH STREET, LONDON, W.1

250, IMPERIAL DRIVE, RAYNERS LANE, HARROW



UNIVERSITY OF LONDON

**INSTITUTE OF OPHTHALMOLOGY, JUDD STREET, LONDON, W.C.1**Associated with  
**MOORFIELDS EYE HOSPITAL**

Qualified Medical Practitioners may enter the practice of the Moorfields Eye Hospital at any time, and are, on certain conditions, eligible for appointment as Clinical Assistants.

**Systematic courses suitable for Diploma and Degree examinations begin in OCTOBER and MARCH of each year.***For further particulars apply to the Dean.***EXAMINATION OF PATHOLOGICAL OPHTHALMIC SPECIMENS.**

The Department of Pathology of this Institute will examine and report on certain specimens sent by Ophthalmic Surgeons. Specimens sent from Hospitals for further examination should be sent through the particular Hospital Pathologist.



For particulars apply to the Director of the Department of Pathology.

**ROYAL EYE HOSPITAL****ST. GEORGE'S CIRCUS, S.E.1**

Telephone: WATERLOO 4477

The clinical practice of the Hospital is open to postgraduate students. There are over 1,450 admissions annually and large daily out-patient clinics.

Courses of advanced tuition are held for the hospital medical staff; other suitable students are admitted as numbers permit.

**The *Wrayflex* in  
OPHTHALMOLOGY****for Photomicroscopy  
and Macrography**Illustrations show  
the camera fitted  
to a microscope  
and attached  
to copying  
stand*Write for 16pp  
illustrated Brochure***WRAY (OPTICAL WORKS) LTD., BROMLEY, KENT**

# MONOGRAPH SUPPLEMENTS

TO BRITISH JOURNAL OF OPHTHALMOLOGY

- No. 3 THE NATURE OF THE INTRA-OCULAR FLUIDS** Price 7/6  
*(Sir Francis Laking Prize 1926-27)*  
 By W. STEWART DUKE-ELDER.  
 140 pages, with Bibliography.
- No. 4 THE NATURE OF THE VITREOUS BODY** Price 5/-  
 By W. STEWART DUKE-ELDER. 72 pages.
- No. 6 RADIANT ENERGY** Price 7/6  
 as a Pathogenic and as a Therapeutic Agent  
 in Ophthalmic Disorders  
*(The Gifford Edmonds Prize in Ophthalmology)*  
 By H. B. STALLARD. 128 pages.
- No. 9 ON THE CORRELATION OF SOME SENSORY AND PHYSIOLOGICAL PHENOMENA OF VISION** Price 5/-  
 By W. D. WRIGHT and RAGNAR GRANIT. 80 pages.
- No. 10 OCULAR ONCHOCERCIASIS, INCLUDING AN INVESTIGATION IN THE GOLD COAST** Price 7/6  
 By HAROLD RIDLEY.  
 58 pages, with 27 illustrations and 3 colour plates.
- No. 11 UVEAL SARCOMATA** Price 10/-  
 By J. D. MARTIN-JONES.  
 94 pages, 22 illustrations.
- No. 12 OCULAR VERTICAL DEVIATIONS** Price 12/6  
 By J. RINGLAND ANDERSON.  
 108 pages, 37 illustrations.
- No. 13 AN EXPERIMENTAL AND CLINICAL STUDY OF THE REACTION OF THE ANTERIOR SEGMENT OF THE EYE TO CHEMICAL INJURY, WITH SPECIAL REFERENCE TO CHEMICAL WARFARE AGENTS** Price 21/-  
 By IDA MANN, A. PIRIE, and B. D. PULLINGER.  
 176 pages profusely illustrated in monochrome and colour.
- No. 14 THE INCIDENCE AND CAUSES OF BLINDNESS** Price 10/6  
 An International Survey, 140 pages.  
 By ARNOLD SORSBY.

Nos. 1, 2, 5, 7 and 8 are out of print.

Publishing Manager


## BRITISH MEDICAL ASSOCIATION

B.M.A. House, Tavistock Square, London, W.C.1



## MAN OF THE MOMENT

Our model is not wearing spectacles. If it were necessary he would undoubtedly be attracted by a Spiller "Statesman". This truly impressive frame is a cunning combination of plasticimshell and rolled gold. It is an exceptionally robust model yet one that is surprisingly light in wear.



The Spiller "Statesman"

# GEORGE SPILLER,

Ltd.

41, WIGMORE ST., W.1

Tele.: WEL. 7471-2

Also at: 57, New Kings Road, Fulham, S.W.6.

THE FIRST OPTICIAN  
IN WIGMORE STREET.  
ESTABLISHED 1873

Tele: REN 6995