


JOHN WEISS & SON, LTD.


17 WIGMORE STREET,
LONDON, W.1

Telephone: MUSEum 7733


INSTRUMENTS FOR CORNEAL SUTURING



B951b




B874



B1053


- B951b.** Suture Forceps, corneal, St. Martin's, *stainless steel.*
- B874.** Forceps, corneal, Barraquer's *do.*
- B1053.** Scissors, miniature, for removal of corneal sutures, etc.



B607



B608



B609

- B607.** Needles corneal, Jameson Evans', 8, 10 and 14 mm.
- B608.** Ditto *do.* Kalt's, 12, 14 and 18 mm.
- B609.** Ditto *do.* Vogt's, 7, 10 and 14 mm.

You are cordially invited to visit our new Showrooms at 17, Wigmore Street, W.1, and inspect a comprehensive range of Ophthalmic Instruments and Apparatus.

CLIFFORD BROWN L^D.

45 WIGMORE STREET, LONDON, W.1

MANUFACTURING DISPENSING OPTICIANS

take pleasure in announcing the opening of a

NEW BRANCH

AT


36 HEATH ROAD

TWICKENHAM

Telephone: POPESGROVE 9543


BRANCHES AT

**EAST SHEEN, EXETER, GLASGOW, KINGSTON, MANOR PARK,
MORDEN, ORPINGTON, PERTH, PUTNEY, REDHILL, REIGATE,
RICHMOND, S. WIMBLEDON, SUNDERLAND, VICTORIA,
WANDSWORTH, WATERLOO, WELLING, WIMBLEDON, WORTHING,
WEST WORTHING**



Hamblin Universal Slit Lamp


This low-priced Slit Lamp, simple and robust, ensures freedom from complex maintenance worries. The microscope and light source, being independent of one another, permit an infinite variety of adjustments. Settings to the chin and head rest are quick and easy. Vertical displacement of the light spot assists in the detection of aqueous flare. As the instrument requires no special table, floor space need not be occupied permanently. The compactness of this Slit Lamp makes it readily transportable.




Illumination is by a 12-volt, 24-watt bulb off a mains transformer (in remote places the Slit Lamp can be run from a car battery). Objectives ($\times 2$ and $\times 4$) and eyepieces ($\times 5$ and $\times 10$) provide a magnification range of 10, 20 and 40 diameters. An adjustable table is available.

THEODORE HAMBLIN LTD

15 Wigmore Street, London, W.1.



- | | |
|--------------------|---------------------|
| 1. BANBURY | 14. LONDON |
| 2. BATH | 15. LYME REGIS |
| 3. BICESTER | 16. MITCHAM |
| 4. BLANDFORD | 17. OXFORD |
| 5. BOSCOMBE | 18. POOLE |
| 6. BOURNEMOUTH | 19. RINGWOOD |
| 7. BRACKLEY | 20. SEVENOAKS |
| 8. BRIDPORT | 21. SHERBORNE |
| 9. CHARD | 22. TONBRIDGE |
| 10. DIDCOT | 23. TUNBRIDGE WELLS |
| 11. DORCHESTER | 24. WEYMOUTH |
| 12. EAST GRINSTEAD | 25. WIMBORNE |
| 13. FARINGDON | 26. YEOVIL |



MELSON WINGATE LIMITED

“We did not want to leave Wigmore Street. We thought it was a nasty blow after so many years. How glad we are, however, that we were forced to do so. Larger, lighter, and all on one floor, our new premises have greatly pleased all who have entered therein.”

Melson Wingate Ltd

DISPENSING OPTICIANS

**18 NEW CAVENDISH STREET
LONDON W.1.**

TELEPHONE: WELbeck 8710



Toledo by Clement Clarke

**Accurate dispensing . . .
individual styling**

Service to the Consultant and his patient is
available through forty-six branches.

CLEMENT CLARKE LTD.

WIGMORE STREET, LONDON, W.1. Telephone: LAngham 2242




Harcourt by Clement Clarke

**Accurate dispensing . . .
individual styling**

Service to the Consultant and his patient is
available through forty-six branches

CLEMENT CLARKE LTD.
WIGMORE STREET, LONDON, W.1. Telephone: LAngham 2242

The Specialist



Introduced in 1955, this modern, well-designed ophthalmoscope is still a first favourite with the specialist, and symbolic of the many first-class instruments of the Keeler range, all readily demonstrated at our West End Instrument Showrooms.



Keeler

INSTRUMENTS THAT


LEAD THE WORLD

39 WIGMORE STREET, LONDON, W.1, ENGLAND

5536 BALTIMORE AVE., PHILADELPHIA, Pa., U.S.A.

the
**CARL
 ZEISS**
 slit
 lamp

is now
 available
 for
 immediate
 delivery
 from stock



RENOWNED FOR ITS BRILLIANT ILLUMINATION,

the Carl Zeiss Slit Lamp has
 the following additional advantages:

- * Rapid change between five different magnifications
- * Great depth of focus
- * Particularly sturdy design gives years of hard service
- * Easily equipped with full range of accessories:
 - [1] Camera equipment. [2] Hruby lens
 - [3] Tutor eyepiece. [4] Gonioscopy equipment
 - [5] Aplanation Tonometer

Illustrated booklet on request from

Degenhardt
 AND COMPANY LIMITED

6 Cavendish Square, London W1 Telephone LANgham 6097


C. W. DIXEY & SON LTD.

19 WIGMORE STREET · LONDON W.1.

ESTABLISHED 1777

TELEPHONE: LANGHAM 1713


**Dispensing Opticians and Manufacturers of
OPHTHALMIC INSTRUMENTS AND APPARATUS**




D 1172



D 1130



D 1205



D 1160

- D 1172 Corneal Graft Knives, Desmarre's.
- D 1130 Corneal Graft Trephines, Franceschetti's.
- D 1205 Suture forceps, Rycroft's, for micro silk.
- D 1160 Corneal Scissors, Castroviejo's.

*Delivery of the above or other instruments from our comprehensive range
can usually be made from stock.*


Catalogue sent on request.



BY APPOINTMENT OPTICIANS TO HER LATE MAJESTY QUEEN MARY

C. W. DIXEY & SON LTD.
19 WIGMORE STREET : LONDON W.1.

TELEPHONE: LANGHAM 1713



INTRODUCING
THE NEW
ULTRA VIOLET RAY
TONOMETER
STERILIZER

PRICE £21
(For use on A.C. mains only)


Sterilization is achieved by placing the barrel of the Tonometer in one of the larger orifices at the top of the sterilizer. The tonometer is thus supported whilst its base and plunger end is irradiated for approximately 40 seconds. Extensive tests have shown that this treatment destroys bacteria and viruses, including *Ps. pyocyanea*. Provision is made for treating two tonometers together, either assembled or with their plungers removed and placed in the smaller orifices.

(This sterilizer was suggested by S. S. Sumner Esq., F.R.C.S., of Preston and is made exclusively for Dixey by Hanovia Ltd.)

Hamblin Cabinet Test Type

.....

No. 804



Reflections from the glass front are eliminated in the redesigned Hamblin Cabinet Test Type. The lighting is now housed internally; it provides the types with good, even illumination. A selection of up to eight Hamblin distance types is mounted on rollers for remote control by cord. Maddox 'Spot' and Duochrome Test are incorporated.

Types with reversed lettering can be fitted for use with a mirror at three metres.

THEODORE HAMBLIN LTD

15 Wigmore Street, London, W.1.

OCULAR PROSTHETICS

MR. STEER WARDMAN

Artificial Eye Specialist

of

Brooks & Wardman

Dispensing Opticians

Artificial Eye Makers

Mr. Steer Wardman will be pleased to accept your referred patients for the fitting of Artificial Eyes and Cosmetic Contact lenses by appointment at:—


24 REGENT STREET — NOTTINGHAM
Telephone: 46309 (2 lines)

91 CORNWALL STREET — BIRMINGHAM
Telephone: Central 0583


28 WATSON ROAD — WORKSOP
Telephone: 3204

J. W. WARDMAN LTD.

54 DUKE STREET — DARLINGTON
Telephone: 66089



My work covers all recent advances in ocular prostheses, magnetic artificial eyes, and facial prostheses.



**DISPENSING OPTICIANS
AND INSTRUMENT MAKERS**

manufacturers and suppliers of

TRIAL LENSES TO BRITISH STANDARDS 3162/1959

PERIMETRY AND SCOTOMETRY EQUIPMENT

ANTERIOR CHAMBER IMPLANTS

REFRACTION ACCESSORIES

CONSULTING ROOM FURNITURE · COLOUR TESTS

ALL-BURIED SOCKET IMPLANTS

AIDS FOR PATIENTS WITH LOW VISUAL ACUITY

DIATHERMY EQUIPMENT · ERISIPHAKE APPARATUS

OPERATING LAMPS · ORTHOPTIC EQUIPMENT

DISTANCE AND NEAR SUBJECTIVE TESTS

OPHTHALMIC CAMERAS


Branches at:

LONDON · BEXHILL · BIRKENHEAD · BOLTON · BRADFORD · BRIGHTON · BRIXTON · CHESTER · COVENTRY
CREWE · DONCASTER · EASTBOURNE · ELLESMERE PORT · HALIFAX · HARROGATE · HASTINGS · HERFORD
HESWALL · HOVE · HUDDERSFIELD · HULL · LEAMINGTON · LEEDS · LEICESTER · LIVERPOOL · LUDLOW
NORTHWICH · NUNEATON · OLDHAM · OSWESTRY · RHYL · RUGBY · RUGELEY · SHREWSBURY
SOUTHPORT · STAFFORD · ST. HELENS · STOKE · STREATHAM · WALLASEY · W. HARTLEPOOL · WORCESTER

INSTRUMENT DEPARTMENT

160 NEW BOND STREET, LONDON, W.1

MAYfair 0128



Hamblin Purvis Streak Retinoscope

. for the rapid and accurate assessment of refractive errors. The brilliant streak of light projected simplifies the location of the axis of astigmatism; the streak clearly shows on the trial frame, whence the angle may be read. Controls for brightness, width and angle of the streak are within easy reach of the fingers.

Complete in case.



THEODORE HAMBLIN LTD


15 Wigmore Street, London, W.1.

HANDINESS is a general characteristic of the Gambs Slitlamp.

Working controls have been grouped around the principal joy-stick control in order to permit one-handed manipulation with no more than flexion of the wrist. The other hand of the operator is left free to hold open the patient's eye or manipulate the fixation device. This arrangement is particularly useful when applying gonioscopy or aplanation tonometry.

Moreover, around the joy-stick, the base has been designed with great care so that an easy positioning of the hand with firm support of the wrist, together ensure very accurate movements. Additional magnification is available by movement of a lever, with no more than a momentary pause in observation. The microscope remains in focus throughout.

Concentric mounting of the lighting and microscope arms also ensures that both essentials remain in focus. Passing of the arms does not involve interruption of observation.



Since 1927 specialized experience in the design and production of ophthalmic instruments has impressed upon us the value of such attention to detail.



GAMBS
140, RUE MAZENOD
LYON
FRANCE

Agents for U. K. :

STERCKS - MARTIN Ltd - 36-56 New Cavendish Street - LONDON W1
SLIT LAMP, complete with table and transf. refer. 850 T £ 310

Please write for fully descriptive literature of our products
SLIT LAMP — OPHTHALMOMETER — PERIMETER
TRANSLATING TABLES — ADJUSTABLE STOOLS

THE FISON BINOCULAR INDIRECT OPHTHALMOSCOPE



This unique instrument is universally preferred for detachment surgery, diagnosis and teaching.

The Reflecting Plate (inset) enables two observers to view a fundus simultaneously, transmitting about 50% of the light incident on it from the patient's fundus and reflecting 50% towards the second observer. This accessory is thus excellent for teaching fundus appearance to students. Being portable, the "FISON" entirely supersedes other demonstrating ophthalmoscopes since it may be carried round wards and clinics.

Expert staff are pleased to demonstrate this and other instruments of the Keeler range at our West End Instrument Showrooms.




Keeler INSTRUMENTS THAT

LEAD THE WORLD

39 WIGMORE STREET, LONDON, W.1, ENGLAND

5536 BALTIMORE AVE., PHILADELPHIA, Pa., U.S.A.

Zonulysin in Cataract



▲ Zonulysin instilled to flood anterior and bathe suspensory ligament.

◀ Suspensory ligament resolved giving easy withdrawal of crystalline lens.

These illustrations are taken from a film demonstrating the use of Zonulysin in cataract. Copies of the film are available on request for loan to Hospitals.

Zonulysin is a highly purified, sterile freeze-dried Alpha-Chymotrypsin, which is particularly indicated for resolving the ciliary zonule.

Maw's Zonulysin is conveniently presented in boxes of 6 ampoules, each containing 0.3 mgm. Alpha-Chymotrypsin together with 6 ampoules each containing 1.0 ml. of Saline diluent.

Zonulysin in the surgery of Cataract has these advantages in a considerable variety of cases.

Reference 1: Boyd, Benjamin, F. (1959-60) *Highlights of Ophthalmology* 3.4. Pages 278-287.

- No retardation of wound healing
- No allergic reaction in thousands of cases
- Zonular fibres fragmented by the enzyme
- Increased facility of extraction
- Decreased number of ruptured capsules

S. MAW SON AND SONS LIMITED

ALDERSGATE HOUSE

BARNET

HERTFORDSHIRE

*The Highest Standards of
Optical Dispensing . . .*

ALFRED

HAWES


and Son Ltd.

120 YEARS OF UNBROKEN SERVICE

The firm of Alfred Hawes was established in 1840, and is actively controlled today by descendants of the Founder.

Up-to-date techniques and this long experience enable us to offer the widest possible service, but we have never departed from the tradition of personal attention to individual requirements.

That tradition is faithfully maintained in all our branches. We look forward to dispensing the requirements of your patients.



HEAD OFFICE:

79-79a, LEADENHALL STREET, LONDON, E.C.3

TELEPHONE: AVENUE 4538/9

WEST END OFFICE: 25, NEW CAVENDISH STREET, LONDON, W.1. Telephone: WELbeck 7862

BRANCHES:

BOW:	<i>Advance</i> 2103	GUILDFORD:	<i>Guildford</i> 4911	ST. PAUL'S CRAY:	<i>Foots Cray</i> 4803
CHISLEHURST:	<i>Imperial</i> 4433	ILFORD:	<i>Valentine</i> 5825	STRATFORD:	<i>Maryland</i> 3393
CROYDON:	<i>Croydon</i> 5779	LEATHERHEAD:	<i>Leatherh'd</i> 3157	WALTHAMSTOW:	<i>Coppermill</i> 5206
DORKING:	<i>Dorking</i> 3187	LEWISHAM:	<i>Lee Green</i> 2245	WOKING:	<i>Woking</i> 944
ELTHAM:	<i>Eltham</i> 4325	LEYTON:	<i>Leytonstone</i> 2109	WALTON-ON-THAMES:	<i>Walton</i> 1832
GODALMING:	<i>Godalming</i> 2088	ROMFORD:	<i>Romford</i> 42208		



SOFT-LITE No. 1 is back again!

THE WORLD'S PREMIER NEUTRAL ABSORPTION LENS

***NOW MADE IN
ENGLAND***

SINGLE VISION LENSES

•
SOLID BIFOCALS

•
FUSED BIFOCALS

•
**UNIVIS BIFOCALS
and TRIFOCALS**

•
**UNIFORT TOUGHENED
LENSES**


•
**SALVOC LAMINATED
SAFETY LENSES**

•

SOFT-LITE 1 at greatly reduced prices, in a wide range of styles, is again freely available. SOFT-LITE 1 is being made in this country to the exact original specification, and a range of single vision lenses in this glass is available from stock. Bifocals and other special lenses are made to order.

Technical specification required by BS 2738, 1956 available on request.

***REDUCED PRICES —
— QUICKER SERVICE***



**U.K. OPTICAL BAUSCH & LOMB LTD.
BITTACY HILL . MILL HILL . LONDON . N.W.7**

BORALINE

decongestive
mildly
antiseptic

ideal in
conjunctival
irritation




sole distributors:

COATES & COOPER LTD. WEST DRAYTON, MIDDLESEX

ZONULOLYSIS

*Eliminating
accident risk
in
Cataract
Surgery*



CHYMAR-ZON

brand of α -chymotrypsin

PRESENTATION

each carton contains:

1 vial

100 Armour Units
lyophilised
 α -chymotrypsin.

1 vial

1 ml. Ringer's
Solution
as diluent.

1 vial

1 ml. Ringer's
Solution for
irrigation.

The inherent complications of extracapsular extraction and the accidents which sometimes occur during or after intracapsular extraction are eliminated by Barraquer's technique.

This technique employs the zonulolytic action of α -chymotrypsin, effectively replacing the mechanical procedures of zonulotomy.

Literature on request



ARMOUR PHARMACEUTICAL COMPANY LTD
EASTBOURNE · SUSSEX · ENGLAND



LINWOOD STRONG
LIMITED

dispensing opticians

SURREY


20 High Street,
Sutton.
Tel: VIGilant 2798.

45 Upper High St.,
Epsom.
Tel: EPSom 2300.

SUSSEX

Chapel Road,
Worthing.
Tel: Worthing 5193.

57b London Road,
Brighton.
Tel: Brighton 29425.



The LAW Distance Type


(Designed at the suggestion of
Mr. Frank Law)

Patients almost invariably read a distance type throughout, starting at the top and ending at the limits of their visual acuity. In the Law Distance Type, the sequence of standards is reversed, beginning with 6/3 and ending with 6/60. A patient is thus impelled to call the smallest letters he can see at the outset, so saving the ophthalmologist's time.

THEODORE HAMBLIN LTD

15 Wigmore Street, London, W.1.

Fashion Frames



* *

*A further example of the spectacle-maker's craft,
combined with the jeweller's. Hand-made frames,
trimmed hand-set rhinestones and pearls.*

Earrings to match.

* *

Made and Designed by

A. E. NEWBOLD, LTD


11b, Hay Hill, Berkeley Square, London, W.1

HYDE PARK 0946/6

**Branch: 1 Broadway Buildings
BROMLEY, Kent**

**Branch: 64 High Street
BECKENHAM, Kent**

★ THE
GAMBS SLIT LAMP
PHOTOGRAPHIC ATTACHMENT



A simple and inexpensive attachment enables the microscope to be replaced by any suitable single lens reflex camera of the interchangeable lens type in a matter of a few moments

THE EXCEPTIONAL ILLUMINATION IN THE GAMBS SLIT LAMP MAKES POSSIBLE **DIRECT COLOUR PICTURES** OF THE ANTERIOR SEGMENT

Thirty-three years of production experience and continuous research contribute to the success of the GAMBS SLIT LAMP of today.


**SLIT LAMP
OPHTHALMOMETER
PERIMETER
TRANSLATING TABLES
ADJUSTABLE STOOLS
ETC.**

**SOLE AGENTS
U. K. & EIRE**

Sterckx Martin
LTD

INSTRUMENT DEPT. 36 NEW CAVENDISH STREET, W.1 Tel: HUNter 0644

New HAAG-STREIT '900' SLIT LAMP & TONOMETER



Aplanation tonometer
showing pressure arm
and prism.

Prism in contact
with cornea.

*The
HAAG-STREIT
'900' Slit Lamp
is the
finest
available
today*

MAIN FEATURES

- Precentred and prefocused bulb, burning cap uppermost.
- Gives 30% increased illumination.
- Slit rotates through 180°. No deviating prism used.
- Unique inclining device provides horizontal *sections*.
- Hruby Lens pre-set and linked to microscope.
- Non-dazzling fixation target focused at far point.
- Aplanation Tonometer attached and removed at will.
- New Fundus Lens for high myopes available.
- Measuring device available.




Please write for literature and full details to:

SOLE AGENTS FOR THE UNITED KINGDOM

CLEMENT CLARKE LTD.

Instrument Division Wigmore St. London W.1



*This is
the
NEW
BRITISH*


Curpax SLIT LAMP

BASIC SLIT LAMP & MICROSCOPE UNIT £150

The portable
model illustrated here
was designed
by us for
Professor Ida Mann
to aid her
Trachoma research.

It is available
in a number of
other combinations.


Made by the
Instrument Division of **CURRY & PAXTON LTD**



**The Stereoscopic
Microscope
is readily
detachable for
RECUMBENT
GONIOSCOPY**

The BASIC UNIT
may be fitted

- To a counterpoise arm
- On Refractor Stand
- On wall bracket
- On a table model
- In a portable box
- A U.V. filter is available
for contact lens work



Patents pending

**Really BRILLIANT SLIT ILLUMINATION
and CONTRAST**

May be used with the patient seated in Ophthalmic chair without chinrest. There is no obstruction of the patient and he cannot draw away from the beam.

CURRY & PAXTON, LTD

Dispensing Opticians, Manufacturers of Ophthalmic Instruments

22, WIGMORE STREET, LONDON, W.1.

**195-199, GREAT PORTLAND STREET,
LONDON, W.1.**

Telephone: LANGham 0123

**INSTRUMENT DEPARTMENT,
120-126, ALBERT STREET, PARKWAY, N.W.1**

Telephone: GULLiver 7939

**PRINCIPAL BRANCHES — BATH, BRISTOL, CHIPPENHAM, CHISWICK, DERBY, GLASGOW,
LINCOLN, LIVERPOOL, LOWESTOFT, NEW CROSS, NEWARK, NORWICH, NOTTINGHAM, RICHMOND,
SCUNTHORPE, SHEPHERDS BUSH, SOUTHPORT, WARRINGTON, WESTON-SUPER-MARE, WOOLWICH.**

Faithfully yours

PANKHURST & BARRETT LTD.

Dispensing Opticians

27, NEW CAVENDISH STREET, LONDON, W.1

250, IMPERIAL DRIVE, RAYNERS LANE, HARROW

GEORGE ODELL LTD.
DISPENSING OPTICIANS

A. H. ODELL, *Managing Director.*
M. ODELL, *Secretary.*

*Members of the Guild of
British Dispensing Opticians*

Makers of spectacles of a fine precision and of a standard worthy of the best work of British craftsmen—a very important but fast-dwindling section of the optical industry.

Frames are made and mounted to the individual measurements of each patient and both frame and lenses periodically checked during the course of production.

The finished spectacles are finally checked by a qualified dispensing optician and whenever possible, personally fitted on every patient.

**31, KING STREET,
ST. JAMES'S, S.W.1**
Whitehall 9987

**36, WELBECK STREET
W.1**
Welbeck 5320

An invitation to subscribe
to this Important New Pergamon International Research Journal

VISION RESEARCH

VISION RESEARCH is an important new international journal devoted to all aspects of seeing. It provides a common meeting ground for research workers in problems of vision, whose separate approaches are primarily biophysical, physiological, biochemical, neurological, psychological or ophthalmological. The emphasis is on human and animal vision rather than on optics as a physical science, though overlapping studies—for example biophysical or physical-chemical—are considered of the highest relevance. Also included are papers in those areas of vision research which are somewhat neglected in the current scientific literature, such as visual perception, experimental ophthalmology, the sensory neurology of vision, the embryology of vision, comparative sensory physiology, physiological optics, the history of visual science etc.

CONTENTS OF VOLUME 1 NUMBERS 1/2—JUST PUBLISHED

D. FARNSWORTH : Let's look at those isochromatic lines again; F. H. G. PITT : Comments on Dean Farnsworth's paper "Let's look at those isochromatic lines again"; W. D. WRIGHT : Comments written to Dean Farnsworth in 1958; D. L. MacADAM : A non-linear hypothesis for chromatic adaptation; H. C. SPERLING : Prediction of relative luminous efficiency from fundamental sensation curves; H. R. BLACKWELL and O. M. BLACKWELL : Rod and cone receptor mechanisms in typical and atypical congenital achromatopsia; F. H. THOMAS, F. L. DIMMICK and S. M. LURIA : A study of binocular colour mixture; L. L. SLOAN : Area and luminance of test object as variables in examination of the visual field by projection perimeters; F. CRESCITELLI : The electroretinogram of the antelope ground squirrel; KATHERINE TANSLEY, R. M. COPENHAVER and R. D. GUNKEL : Spectral sensitivity curves of diurnal squirrels; E. AUERBACH, A. J. BELLER, H. E. HENKES and G. GOLDHAVER : Electric potentials of retina and cortex of cats evoked by monocular and binocular photic stimulation; R. A. WEALE : Notes on the photometric significance of the human crystalline lens; H. KALMUS : The attention of optomotor responses in white-eyed mutants of *musca domestica* and *coelopa frigida*.

Published Quarterly

BOARD OF EDITORS

- DR. T. SHIPLEY, *Department of Ophthalmology, University of Miami College of Medicine, Miami, 36, Florida*
- F. CRESCITELLI, *Department of Physiology, University of California, Los Angeles*
- H. J. A. DARTNALL, *Institute of Ophthalmology, Judd St., London*
- PROF. Y. LEGRAND, *Museum National d'Histoire Naturelle, Laboratoire de Physique, Paris*
- PROF. DR. H. SCHOBER, *Institut für medizinische Optik, der Universität, München*
- PROF. A. SORSBY, *Department of Ophthalmological Research, Royal College of Surgeons of England, London*

Assisted by an International Honorary Editorial Advisory Board

Manuscripts for editorial consideration should be sent to the Member of the Board of Editors most conveniently located.

A copy of this journal and details of subscription rates gladly sent on application



PERGAMON PRESS

OXFORD LONDON NEW YORK

Headington Hill Hall, Oxford.

4 & 5 Fitzroy Square, London W.1.

122 East 55th Street, New York 22, N.Y.

UNIVERSITY OF LONDON

INSTITUTE OF OPHTHALMOLOGY, JUDD STREET, LONDON, W.C.1

Associated with
MOORFIELDS EYE HOSPITAL

Qualified Medical Practitioners may enter the practice of the Moorfields Eye Hospital at any time, and are, on certain conditions, eligible for appointment as Clinical Assistants.

Systematic courses suitable for Diploma and Degree examinations begin in OCTOBER and MARCH of each year.

For further particulars apply to the Dean.

EXAMINATION OF PATHOLOGICAL OPHTHALMIC SPECIMENS.

The Department of Pathology of this Institute will examine and report on certain specimens sent by Ophthalmic Surgeons. Specimens sent from Hospitals for further examination should be sent through the particular Hospital Pathologist.

For particulars apply to the Director of the Department of Pathology.

ROYAL EYE HOSPITAL**ST. GEORGE'S CIRCUS, S.E.1**

Telephone: WATERLOO 4477

The clinical practice of the Hospital is open to postgraduate students. There are over 1,450 admissions annually and large daily out-patient clinics.

Courses of advanced tuition are held for the hospital medical staff; other suitable students are admitted as numbers permit.

UNIVERSITY OF SYDNEY**CHAIR OF EYE HEALTH**

Applications are invited for the recently endowed Chair of Eye Health.

Salary will be at the rate of £A4,930 per annum, plus cost of living adjustments (at present £A42 p.a.). There is retirement provision under either the Sydney University Professorial Superannuation Scheme or the New South Wales State Superannuation Scheme.

Under the Staff Members' Housing Scheme, in cases approved by the University and its Bankers, married men may be assisted by loans to purchase a house.

The Senate reserves the right to fill the Chair by invitation.


A statement of conditions of appointment and information for candidates may be obtained from the Secretary, Association of Universities of the British Commonwealth, 36 Gordon Square, London, W.C.1.

Applications close, in Australia and London, on 15th August, 1961.

AN AID TO LOW VISUAL ACUITY

WRAY PANORA

LONDON



Can be supplied under the National Health Service if prescribed as clinically necessary by an ophthalmologist in the National Eye Service.

Magnification $\times 8$.

Object Glass 21 mm.

Folds flat into soft leather case

$2\frac{3}{4}'' \times 1\frac{7}{8}'' \times \frac{3}{4}''$.

Details from

WRAY (OPTICAL WORKS) LTD

BROMLEY, KENT

Telephone RAV 0112

W. CALLAGHAN & S. P. HOLLOWAY LTD.

ACCURATE & PROMPT DISPENSING

PERSONAL ATTENTION BY
THE PRINCIPALS

14A, HENRIETTA PLACE

(Corner of Welbeck Street)

LONDON · W.1

Telephone Welbeck 2703

Selection


From some of the
World's Finest Spectacleware



A pretty and practical new German style.



... dispensed
with the inherited
skill of 86 years'
experience ...



... plus
attention
that is really
personal ...

... at

SPILLERS

GEORGE SPILLER,
Ltd.

41, WIGMORE ST., W.1

Tele.: WEL. 7471-2

Also at: 57, New Kings Road, Fulham, S.W.6. Tele: REN 6995
457A, Finchley Road, N.W.3. Tele: HAM 5524

THE FIRST OPTICIAN
IN WIGMORE STREET.
ESTABLISHED 1873