Sno-Pilo

OPPORT

pilocarpine you can always be sure of.



Sno-Pilo is presented in a specially designed, easy-to-use, plastic dropper bottle containing 10ml pilocarpine hydrochloride eye drops. There are five strengths available: 0.5% w/v, 1.0% w/v, 2.0% w/v, 3.0% w/v and 4.0% w/v. PL 0033/0064 – 0068. Full prescribing information is available. Smith & Nephew Pharmaceuticals Ltd., Bessemer Road, Welwyn Garden City, Hertfordshire, United Kingdom.



The Isopto* vehicle



CONTROLL OF THE PROPERTY OF TH

The ideal ophthalmic delivery system which

Prolongs corneal contact time¹ Produces a greater reduction in IOP²

Lubricates and soothes.

An effective stable product with precisely the strength you specify, every bottle, every batch which gives a flexible range of strengths 0.5%, 1%, 2%, 3%, 4% and 6% together with the benefits of the Alcon Drop-Tainer enabling your patient to instil drops with accuracy and ease.

Linn M.L. and Jones L.T.
 Amer.J. Ophthal. 65. 76. 1968
 Haas J.S. and Merrill D.L.
 Amer. J. Ophthal 54. 21. 1962

Alcon

ALCON LABORATORIES (U.K.) LIMITED Imperial Way Watford Hertfordshire England WD2 4YR Felephone Watford 46133 Telex 923709 Cable Alcon Watford




pilocarpine you can always be sure of.



Sno-Pilo is presented in a specially designed, easy-to-use, plastic dropper bottle containing 10ml pilocarpine hydrochloride eye drops. There are five strengths available: 0.5% w/v, 1.0% w/v, 2.0% w/v, 3.0% w/v and 4.0% w/v. PL 0033/0064 – 0068. Full prescribing information is available. Smith & Nephew Pharmaceuticals Ltd., Bessemer Road, Welwyn Garden City, Hertfordshire, United Kingdom.




It's the Originall It's the Heag-Streit Perimeter 940



and now... delivery from stock



Originals have value in the healing arts as well as in the fine arts. The GOLDMANN PERIMETER 940 by Haag-Streit is a prime example.



A product of the world's ophthalmological instrument leader, the Goldmann Perimeter provides high accuracy and economy of time.

The unit standardizes target size, luminosity of target and background and all other examination variables as well.

With its unique, patient-activated recording device, the Model 940 reduces errors and simplifies target presentation control.

There are two original Goldmann models, the 940-K7 for kinetic perimetry and the 940-ST for static and kinetic perimetry.


If you are ready to make the lifetime investment that owning a true original represents, contact Clement Clarke International Ltd. They will be glad to send full information on the Goldmann Perimeter and our many other fine products.



Clement Clarks International Ltd.

15 Wigmore Street, London WIHSLA

Representative in Eire:



TWENTY-FOUR GOOD REASONS FOR USING MINIMS**

MINIMS-STERILE, SINGLE-DOSE EYE DROPS

There are twenty-four presentations in the Minims single-dose range. Each one selected for its usefulness in day to day clinical situations.

Minims solutions are sterile and, because of the over-wrap and special cap closure, remain so right up to the moment of application. No preservatives are included in the formulations, so Minims may be used in open eye surgery.

Full prescribing information is available from Smith & Nephew Pharmaceuticals Limited Welwyn Garden City, Herts, England



Alcon

Pneumatonograph



THE FIRST tonograph which measures intraocular pressure by applanating the eye ... the Alcon Applanation Pneumatonograph. It also brings a new dimension to the rapidly developing field of ocular pulse studies.

The Alcon Applanation Pneumatonograph encompasses the following unique features:

- ■Soft tip, self-aligning applanating sensor
- Accurate recorder reading directly in mm Hg
- Remotely positionable magnetic digital display
- ■Modular Design . . . allowing for expansion of system
- ■Two or four minute tonography
- ■Heat stylus ... no ink required
 The Alcon Applanation Pneumatonograph is
 ... exciting ... versatile ... and opens a new
 horizon for glaucoma management.

"The Pneumatonograph has the advantage of being reproducible in different postures: the influence of posture on intraocular pressure can therefore be more exactly determined. Applanation pneumatonography is more versatile as it is portable and does not require a slit lamp or fluorescein, and a permanent written record is available."

"M. R. JAIN and V. J. MARMION BRIT. J. OPHTHAL. (1976) 60, 107".

For further information apply:-

Alcon Laboratories (U.K.) Ltd., Ophthalmic Division,

Imperial Way, Watford, Herts WD2 4YR Tel: 46133