



TWO USEFUL SPECTACLE LOUPES



THE BISHOP-HARMAN LOUPE




THE MAYOU LOUPE

The two Loupes illustrated above are cheap and simple forms, giving adequate magnification and leaving a large free field. They can be worn over ordinary spectacles and are the lightest and simplest of all binocular Loupes. The upper illustration depicts the Bishop-Harman Loupe in which the circular lenses are in effect the two halves of a large sphere with the consequent prism base in. The alternative Mayou Loupe below is preferred by some as being more easy to wear and having a slightly larger field laterally.

**THEODORE
HAMBLIN LTD**
DISPENSING OPTICIANS
15 WIGMORE STREET,
LONDON, W.1.

The PANTOSCOPE




39, WIGMORE STREET,
LONDON, W. I.



617 S. 52nd STREET,
PHILADELPHIA 43, Pa.



*Completely
absorbed*



Gelatin Sponge A & H, a sterile absorbable hæmostatic, has application in all branches of surgery, from the simple epistaxis encountered in general practice to the more complex operations of neurosurgery.

It minimizes blood loss when adequate hæmostasis cannot be obtained by ligature, and will help to increase the safety of operative surgery in many fields.

Gelatin Sponge A & H is pliable and may be moulded easily to any shape or size to meet the requirements of the occasion.


Supplied sterile in glass containers in strips 2 cm. \times 6 cm. \times 0.7 cm., in boxes of six tubes each containing one piece; in strips 10 cm. \times 20 cm. \times 0.1 cm. in tubes containing one piece, and thin wafers 2 cm. \times 2 cm. \times 0.1 cm. in boxes of six tubes each containing six pieces.



GELATIN SPONGE A & H

An absorbable hæmostatic

Literature and further particulars will be sent on request.



new advance
in antibiotic
topical therapy

Myciguent Ophthalmic Ointment, a new Upjohn product, contains as its active principle neomycin sulphate, newest of broad-spectrum antibiotics of proved usefulness.

Rapidly effective in the treatment of sties, marginal blepharitis, and infections of the conjunctiva and eyelids, Myciguent Ophthalmic Ointment is active against a wide range of both Gram-positive and Gram-negative organisms.

Neomycin is well tolerated, with a sensitizing index which is exceedingly low even in patients known to be sensitive to other antibiotics.

Myciguent*
 OPTHALMIC OINTMENT

Containing 5 mg. neomycin sulphate per Gm. Supplied in 1-drachm tubes.

Myciguent is also available in ointment form for topical therapy of specific or mixed cutaneous infections.

*TRADE MARK

MYCIGUENT OINTMENT—


Supplied in 1 oz. tubes, 5 mg. per Gm.



FINE PHARMACEUTICALS SINCE 1886

UPJOHN OF ENGLAND LTD., 4 ALDFORD STREET, PARK LANE, LONDON, W.1

GROSVENOR 5561



Artificial Eyes in Glass or Plastic.
Embodying All Modern Improvements,
Specially Made to Surgeons' Instructions and
Patients' Requirements.

Appointments may be made at:—

2 BLOOMSBURY PLACE, LONDON, W.C. 1

(Nearest underground station Holborn)

Telephone: Mus. 0908

OPHTHALMIC SURGICAL INSTRUMENTS

by


STEARMAN

Stearman Optical Company, Ltd., well known for the British Refracting Unit and other sight-testing equipment now introduces ophthalmic surgical knives of the highest quality. We are able to give reasonable delivery on knives, forceps, scissors, etc.



SERVICING DEPT.

Repairs to all kinds of ophthalmic surgical instruments executed. High class workmanship. Enquiries invited.




Price lists gladly sent upon application.

STEARMAN OPTICAL CO. LTD.


10 ST. CROSS STREET, HATTON GARDEN, LONDON, E.C.1

Telephone: CHAncery 5485-6

THE LISTER-MORTON OPHTHALMOSCOPE
WITH BATTERY IN HANDLE AND WITH
ALTERNATIVE HANDLE FOR USE
WITH MAINS



No. 495
with battery
in handle.



No. 0495/9
with handle for
use on mains with
transformer or
resistance.

Introduced in 1928 Hamblin's Lister-Morton Ophthalmoscope remains the Standard amongst electric ophthalmoscopes with a performance that has yet to be surpassed. The lighting system has been designed to cut down corneal reflex to a minimum and to eliminate reflections, thus making it possible to carry out a direct examination of the macular region in a fully lighted room without using mydriatics.

The battery model being self-contained is ideal for bedside examination, but operation from the mains is usually preferred in the consulting room.

A good arrangement for a busy Surgeon is to use the battery model for visiting purposes and to have an interchangeable mains handle complete with lamp and condensing system for consulting room purposes.

**THEODORE
HAMBLIN LTD**
DISPENSING OPTICIANS
15 WIGMORE STREET
LONDON, W.1.

GEORGE SPILLER LTD.

DISPENSING OPTICIANS

41 WIGMORE STREET, W.1.

'Phone: Welbeck 7471-2.



We have maintained that traditional craftsmanship in association with advances in optical dispensing which enables us to offer you an unrivalled optical service

SUPREME SERVICE SINCE 1873

—SPILLER