


# Retinophot 211

The first fundus camera with quick changer for viewing angle and magnification – an improvement in retinal photography.



These are the features of the Retinophot 211

- Viewing angles of 40°, 29° and 14.5° corresponding to 15x, 20x, and 40x magnification.
- For serial photography built in automatic mechanism with 0.5 sec film rate, automatic film transport and simultaneous exposure of patients data. Film rates of 1 sec and 2 sec also possible.
- Photographic scale 1.9x, 2.6x, 5.3x.
- Precentered 6V 20W Halogen lamp and precentered flash lamp for flashes up to 120Ws.
- Compact design and counterbalanced movement of the retinal camera: 60° vertical axis, 10° horizontal axis.
- Accessories as second observer eyepiece, Polaroid and TV adapter.

If you require further information or a personal demonstration of the Retinophot 211, please let us know by completing the coupon below.

**C.Z. Scientific Instruments Ltd., Carl Zeiss Jena Ltd.,**  
P.O. Box 43, 2 Elstree Way, Borehamwood, Herts, WD6 1NH.  
Telephone 01-953 1688.



## Precision Optics from VEB Carl Zeiss JENA



Please send me full details of the Retinophot 211  
Demonstration  
(tick as appropriate)

Name \_\_\_\_\_

Company \_\_\_\_\_

Address \_\_\_\_\_

\_\_\_\_\_

BJO



SMITH & NEPHEW PHARMACEUTICALS LTD. OPHTHALMOLOGICAL PRODUCTS

# Simplene®

Adrenaline B.P.

EFFECTIVELY REDUCES INTRA-OCULAR PRESSURE IN OPEN ANGLE, PRIMARY AND SECONDARY, GLAUCOMA

Simplene, supplied as a sterile ophthalmic solution of adrenaline B.P., is presented in a 7.5ml plastic dropper bottle, designed to ensure patient convenience: available in two strengths, Simplene 0.5% and Simplene 1.0%

#### DOSAGE AND ADMINISTRATION:

Adult: One drop to be instilled into the eye once or twice daily.

Children: At the discretion of the physician.

#### CONTRA-INDICATIONS, WARNINGS, etc.:

Simplene should not be used when the diagnosis of open angle glaucoma has not been verified. It is contra-indicated in patients with a narrow-angle because pupillary dilation may precipitate angle-closure glaucoma.

Occasionally patients may complain of orbital discomfort or red eye. Rarely, headache, irritation and local skin reactions occur. As with other adrenaline preparations, melanosis may occasionally

occur, but this has no pathological significance. Systemic effects are rare but include tachycardia, extrasystoles, and elevation of blood pressure.

When used in conjunction with miotics, Simplene should follow the miotic after an interval of 5-10 minutes.

PL 33/57,72.



SMITH & NEPHEW PHARMACEUTICALS LTD.

WATFORD, HERTS. AL9 1DA, ENGLAND