

Supplemental Table 1 Results for repeatability for the fundus blood flow perfusion examination with OCTA (N=9)

Blood flow perfusion	Mean of difference for two continuous tests (95%CI)	t value^a	P value
Superficial retinal capillary	0.065 (-0.140, 0.270)	1.855	0.101
Deep retinal capillary	0.172 (-0.480, 0.825)	1.554	0.159
Choroidal capillary	-0.029 (-0.658, 0.601)	-0.267	0.796

a. Paired T tests were applied to compare the values for two continuous tests of OCTA.