



OPTICIANS
BY APPOINTMENT TO
HIS MAJESTY THE KING

**THEODORE
HAMBLIN LTD
LONDON, W.1.**

MAKERS OF SPECTACLES TO
SURGEONS' PRESCRIPTIONS ONLY.



OPTICIANS
BY APPOINTMENT TO
HER MAJESTY THE QUEEN

YOUR OLD OPHTHALMOSCOPE

CAN BE CONVERTED INTO A LISTER-MORTON OPHTHALMOSCOPE

The introduction of the Lister-Morton Ophthalmoscope a few years ago completely revolutionized accepted ideas as to the amount and quality of the illumination that could be expected from a self-luminous instrument.

The lamp and system of condensers in the Lister-Morton Ophthalmoscope afford an illumination which hitherto had not been deemed possible.

Surgeons who may possess a good Morton electric ophthalmoscope are often disappointed with the quality of the illumination from their own instrument after having tried a Lister-Morton Ophthalmoscope.

Theodore Hamblin, Ltd., can convert such an instrument into a Lister-Morton type for £6 15s. if it is a Luminous Morton and they can use the existing handle, and for £7 15s. if a new handle has to be supplied.

The resulting instrument is identical in appearance and utility with a new Lister-Morton Ophthalmoscope. The price of a new instrument is £12 12s.

PRICES—

Lister-Morton Ophthalmoscope
£12 12 0

Converting from own Instrument
£6 15 0 or £7 15 0

15. WIGMORE STREET, LONDON, W.1
AND AT
MANCHESTER, LIVERPOOL, SHEFFIELD, LEEDS, EDINBURGH,
NEWCASTLE-UPON-TYNE, BOURNEMOUTH, WINDSOR, KINGS LYNN.

ii.


ADVERTISEMENTS

Curry & Paxton, Ltd.
DISPENSING OPTICIANS

195 - 199, G^T PORTLAND ST^E and
22, WIGMORE ST^E LONDON, W.1

Telephones :
Welbeck 0123
Langham 2873.

Telegrams:
"Optician, Wesdo, London"




THE MADDOX VISU-AXIMETER.

(Prov. Patent 3351 29).

AN Apparatus for the determination of precise lens centering based on the visual lines. Curry & Paxton employ this instrument in all cases of high refractive errors, cases of facial asymmetry, and for the accurate fitting of bifocal glasses.

The New Keeler Equipment for Heating the Maddox Eye Warming Pad.



THIS method of heating an eye pad is perfectly safe as it forms a unit independent of the ordinary high-voltage house supply, and removes any possibility of shock to the patient or operator should a wire become exposed. In addition to the heating pad, there is provision for a cautery needle.

This equipment, together with suitable accumulators for supplying the heat to the eye-warming pad, or cautery needle, is mounted in a strong, well-finished portable teak case.

The capacity of the accumulators is 70 ampere hours and they will remain fully charged for one week.

Sole Manufacturer—

C. DAVIS KEELER,

Ophthalmic Instrument Manufacturer and Oculists' Dispenser,

47, WIGMORE STREET, LONDON, W.1

BRANCHES—PLYMOUTH, TRURO and PENZANCE.

COMPLETE EQUIPMENT FOR EYE MUSCLE TRAINING

The Synoptophore

A precision instrument used for training the fusion sense and correcting manifest squint in young children, and exercising old patients for the correction of a phoria.

The Training Stereoscope

Adapted for applying the new "Eye and Hand" method of training an amblyopic eye as suggested by Dr. E. E. Maddox.

The Variable Prism Stereoscope

Stereoscopic exercises for manifest and latent squint, measuring muscle balance.

The Coloured Light Phoria Demonstration Apparatus, etc., etc.

*Booklet giving full description and brief technique sent on
application to :*

CLEMENT CLARKE
dispensing opticians^{LTD}
16 WIGMORE ST.
LONDON W.1
