

# A brighter start for dry eyes...



## begins last thing at night

Used at night Lacri-lube soothes,  
protects and lubricates dry eyes  
for undisturbed sleep and  
comfort on waking.

# LACRI-LUBE®

petrolatum mineral oil OINTMENT

**Lacri-Lube Prescribing Information.** Presentation: Sterile, bland, non-medicated ophthalmic ointment for topical administration to humans, containing white petrolatum mineral oil, non-ionic lanolin derivatives with chlorbutanol 0.5% as a preservative. **Uses:** Useful as adjunctive therapy to lubricate and protect the eye in conditions characterised by exsive keratitis, decreased corneal sensitivity, recurrent corneal erosions, and keratitis sicca. **Dosage and administration:** For topical administration. Pull lower lid down to form pocket. Apply small amount as needed. **Contra-indications etc:** No known contra-indications. **Pharmaceutical precautions:** Store away from heat. To avoid contamination during use, do not touch tip to any surface. **Legal Category: P.** **Package quantities:** Available in 3.5g ophthalmic tubes. **Basic NHS Cost:** (as at May 1987) £1.76. **Product Licence Number:** PL 0426/0041. **Product Authorisation Number:** PA 148/38/1 (Allergan Pharmaceuticals (I) Limited, Westport, Co Mayo).

## At the end of the day there's nothing quite like it for dry eyes




 **ALLERGAN PHARMACEUTICALS**

**SCIENCE WITH VISION**

Allergan Limited,  
High Wycombe, Bucks. HP12 3SH

Full prescribing information is available on request. ®denotes registered trademark

# TOPCON TRC - NW3 NON-MYDRIATIC RETINAL CAMERA



***In keeping with our philosophy of supplying only top quality equipment, Keeler are now the exclusive distributors for TOPCON.***

The fully integrated design of the Topcon TRC-NW3 from Keeler makes it the perfect choice for top quality retinal photography. Featuring an omnidirectional joystick control and a choice of 45° or 20° coverage at a touch, it is fast and easy to use, requiring minimal operator training. The TRC-NW3 produces exceptionally sharp, large diameter 74mm photographs on fine quality Polaroid film.

- Monitor and controls in line with patient to simplify operation.
- Operation possible without complete darkness, as only 4mm pupil dilation is required.
- Full-frame 74mm fine resolution Polaroid film used.
- Auto-wind Polaroid film back incorporates LED operating display.
- Adjustable 3-step flash system allows choice of intensity.
- Automatic internal fixation lamps shift for left/right coverage.
- Green external fixation available.



Please arrange for your Representative to call	<input type="checkbox"/>
Please send me further information	<input type="checkbox"/>
Please use capitals and address your reply to Laura Darragh	
NAME _____	
ADDRESS _____	
POSTCODE _____	
TELEPHONE _____	FAX NO. _____

## Keeler Ltd


Clewer Hill Road, Windsor, Berkshire SL4 4AA  
Telephone: (0753) 857177 Telex: 847565/846588 Fax: (0753) 857817

**London Centre**  
Keeler Martin and Tomkins  
38 Margaret Street  
Cavendish Square  
London W1A 4NS  
Telephone: 01 935 8512


**Scotland Centre**  
25 Deerdrykes View  
Westfield Estate  
Cumbernauld G68 9HN  
Telephone: (0236) 721214  
Fax: (0236) 721231

# No frills No gimmicks . . .


*just the best  
all round  
value for  
money!*




260U




250U



200U



150U



320U


**RAYNER** . . . the name synonymous with IOLs

**Rayner Intraocular Lenses Limited**

1-2 Sackville Trading Estate, Sackville Road, Hove, E.Sussex BN3 7AN, ENGLAND  
Tel. Brighton (STD 0273) 205401 Telex 877744 (RILSTE G) Fax 0273 24623

# Haag-Streit BQ Slit Lamp


The Haag-Streit BQ Slit Lamp is descended from the thoroughbred 900 Slit Lamp, and introduces a new concept in Slit Lamp microscopy with the Stereo-Variator. This allows much wider stereoscopic viewing of the fundus than was previously possible. The bio-microscope had also been re-designed to allow more flexibility whilst maintaining the exceptional optical quality for which the Haag-Streit range has become renowned.



**Clement Clarke International Ltd**

15 Wigmore Street London W1H 9LA England.  
Tel 01-580 8053 Tlx 298625 Fax 01-409 1062

DIXEY



- D276** Govan's Depressor/Marker, Scleral
- D1228** Birks Micro 'Colibri' Forceps
- D1226/22** Beaked (Colibri) Micro Forceps:

**DIXEY INSTRUMENTS LIMITED**

Ophthalmic Instruments & Apparatus  
19 Wigmore Street  
London W1H 9LA England  
Telephone: 01-580 1713 Telex: 268312

# In a hurry for Healonid®?

## Running low on lenses?

At Pharmacia, our reputation for innovative ophthalmic products is matched only by our commitment to efficient service. To satisfy the increasing demand for urgent deliveries of Healonid® and IOLs, we offer the following:

- Datapost delivery
- Courier services
- Personal delivery
- London sales
- Tel: 01-935 5555
- Optional purchasing schemes
- Training courses
- Meeting information service



Contact Birgitta Keenan, our Customer Services Administrator on **0908 670413** for further details.

**Pharmacia Limited**  
Ophthalmics Division  
Pharmacia House Midsummer Boulevard  
Milton Keynes MK9 3HP  
Tel: 0908 661101 (24 hour answering service)



**Pharmacia**

Ophthalmics